DCA® Systems

HEMOGLOBIN A1C SELF-TEST USE

INTENDED USE

The DCA® Hemoglobin A1c assay (HbA1c) provides a convenient quantitative method for determining hemoglobin A1c in blood. The measurement of hemoglobin A1c concentration is recommended for monitoring the long-term care of persons with diabetes. For self-test use.

PRINCIPLES OF PROCEDURE

The American Diabetes Association (ADA) recommends measurement of hemoglobin A1c levels 2-4 times per year, less frequently in patients with stable control.¹

Studies have shown that regular monitoring of hemoglobin A1c values can lead to changes in diabetes treatment and improved metabolic control. This is indicated by decreased hemoglobin A1c values over time.^{1,2}

The DCA system and reagent cartridges measure the HbA1c values in blood. Once the reagent cartridge containing the blood sample is loaded into the DCA test compartment, the blood and reagents are mixed to activate the reaction. The test time is approximately 6 minutes. At the end of the test, the result is automatically calculated by the DCA System and displayed in mmol/mol values, the International Federation of Clinical Chemistry (IFCC) standard.

MATERIALS PROVIDED


- 10 reagent cartridges
- 10 capillary holders
- Calibration card
- Package insert

ADDITIONAL ITEMS NEEDED

- DCA Vantage® System
- Lint-free tissue
- Gloves
- 10% bleach solution (0.525% sodium hypochlorite)
- Gauze pad or cotton ball
- Bandage
- Optional: Quality Control materials, Siemens Healthcare Diagnostics recommends the use of the DCA Normal and Abnormal Control kit.

REAGENT CARTRIDGE COMPONENTS

FIGURE 1. REAGENT CARTRIDGE

- 1 Pull-tab (pull to release buffer from tray)
- 2 Cartridge removal tab
- 3 Capillary holder
- 4 Optical window

Antibody Latex: HbA1c-specific mouse monoclonal antibody adsorbed onto latex particles. 2.5% w/v antibody-latex in 10 mM glycine buffer; 16% w/v nonreactive ingredients (10 µL dried in each reagent cartridge).

Agglutinator: 0.005% w/v poly (aspartic acid) polymer covalently attached to the HbA1c hapten in 20 mM sodium citrate buffer containing 0.1% w/v bovine serum albumin and 1% w/v nonreactive ingredients (10 µL dried in each reagent cartridge).

Buffer Solution: 8.1% w/v lithium thiocyanate, 0.01% digitonin in 200 mM glycine buffer (0.6 mL in each reagent cartridge).

Oxidant: 1.5% w/v potassium ferricyanide in water with 21% w/v nonreactive ingredients (10 µL dried in each reagent cartridge)

- H411 - Toxic to aquatic life with long lasting effects.
- P273 - Avoid release to the environment.
- P391 - Collect spillage.
- P501 - Dispose of contents and container in accordance with all local, regional, and national regulations.

Safety data sheets (MSDS/SDS) available on www.siemens.com/poc

WARNINGS AND PRECAUTIONS

- For optimal performance follow the routine maintenance and cleaning steps outlined in the operator's guide.
- Don't ingest or use internally in any way.
- Keep reagent cartridge and capillary holder away from children.
- When handling the reagent cartridge, don't touch or otherwise contaminate the optical window as this may generate error codes during testing.
- Reagent cartridges and capillary holders are single use only; never perform a second test using the same reagent cartridge or capillary holder.
- To prevent injury, don't forcibly remove a reagent cartridge from the instrument. Ask a healthcare professional for assistance if the problem can't be solved.
- Discard the lancet, reagent cartridge and capillary holder as directed by a healthcare professional.
- To prevent contamination, use gloves when interacting with the instrument.
- Use gloves when assisting someone else in sample collection.
- If test isn't performed by a healthcare professional, disinfect the front of the instrument before each use. Disinfect using a cloth or wipe dampened with 10% bleach solution (0.525% sodium hypochlorite). Don't spray directly onto the instrument.

STORING AND USE

Temperature Indicator

Upon receipt of this kit, check the temperature indicator located on the front of the carton. If the indicator has turned red, don't use the reagent cartridges.

Note The time and date you received the kit, and for assistance in obtaining a replacement kit, contact your local technical support provider or distributor.

Storage

- 8°C Store reagent cartridges refrigerated at 2–8°C (36–46°F).
- 30°C Capillary holders may be stored at room temperature at 15–30°C (59–86°F) or refrigerated at 2–8°C (36–46°F).
- Reagent cartridges can be kept for up to 3 months at room temperature any time before the EXP expiration date.
- Record the date the carton was placed at room temperature on the carton.

Stability of Reagent Cartridges

- After opening the foil pouch, the reagent cartridge must be used within 60 minutes.
- Don't use reagent cartridges after the last day of the expiration month.


GETTING READY FOR TESTING

1. Allow the reagent cartridge to warm up to room temperature if refrigerated.
 - 10 minutes if in an unopened foil pouch.
 - 5 minutes if removed from the foil pouch.
2. Open the foil pouch and tear down from the corner notch until the entire long side of the pouch is open (See Figure 2). Don't use scissors to open the foil pouch.

Important The reagent cartridge must be used within 60 minutes after opening the foil pouch.

Important Discard the reagent cartridge if it's damaged, the pull-tab is loose or missing, the desiccant is missing, or if loose desiccant particles are found inside the foil pouch.

FIGURE 2. OPEN REAGENT CARTRIDGE

3. Remove the capillary holder by peeling back the white plastic film from the clear plastic blister. Don't push the capillary holder through the plastic (See Figure 3).


FIGURE 3. REMOVE CAPILLARY HOLDER

Discard the plastic capillary holder if either of the following is missing:

CAPILLARY HOLDER

- 1 Glass capillary
- 2 Absorbent pad

TEST PROCEDURE

STEP 1: COLLECTING A BLOOD SAMPLE

1. Use a single-use lancing device to collect a blood sample.
 - Clean the finger with an alcohol wipe, or wash hands in warm, soapy water. Dry the fingertip.
 - Prick the finger using a lancing device according to the manufacturer's instructions.
 - Gently squeeze from the base of the finger to form a round drop of blood. If the blood smears or runs, wipe it off with a tissue and gently squeeze another round drop of blood.
2. Fill the capillary holder using blood from a finger stick.

Important Analysis must begin within 5 minutes of filling the capillary holder with the sample.

- Hold the capillary holder at an angle.
- Touch only the tip of the capillary to the small drop of blood on the finger until the capillary fills.

STEP 2: CLEANING THE CAPILLARY

1. Wipe the side of the capillary with a lint-free tissue.
2. Check the glass capillary for bubbles.
 - If bubbles are obvious, discard the capillary holder and repeat the procedure using a new capillary holder.

STEP 3: INSERTING THE CAPILLARY HOLDER

1. Insert the capillary holder into a reagent cartridge with the flat side towards the cartridge, until the holder snaps into place.

STEP 4: SCANNING THE REAGENT CARTRIDGE

1. Hold the reagent cartridge so that the barcode faces to the right.
2. Insert the reagent cartridge into the instrument's track above the black dot.
3. Quickly and smoothly, slide the reagent cartridge down the track. A beep and a screen change indicate a successful scan. Repeat steps for scanning cartridge if this wasn't successful.

STEP 5: INSERTING THE REAGENT CARTRIDGE

1. Open the compartment door and insert the reagent cartridge, with the barcode facing to the right, until you hear a click.
2. Hold the cartridge in place, and using a smooth, slow, continuous motion, pull the flexible pull-tab completely out of the reagent cartridge.
3. Close the compartment door and dispose of the flexible pull-tab. Five seconds after the door is closed, a beep sounds and the test begins.
4. Follow any on-screen prompts.

STEP 6: RECORDING THE RESULTS

1. Wait until the test is complete (approximately 6 minutes). Don't attempt to open the compartment door or remove the reagent cartridge while the test is in progress.
2. Print or record the results from the screen.

STEP 7: REMOVING THE REAGENT CARTRIDGE

1. After viewing and recording the results, select **Next** and open the compartment door.
2. Push and hold down the beige button on the right side of the compartment.

- With your left hand, gently push the clear tab on the reagent cartridge to the right to unlock.
- Pull reagent cartridge out of the compartment.
- Discard the reagent cartridge to prevent injury or possible contamination.

WHAT TO DO WITH THE RESULTS

Note Don't make any decision of medical relevance without first consulting your healthcare provider.

Measuring Range 4–130 mmol/mol

If results are less than 4 mmol/mol, the instrument will display results as < 4 mmol/mol.

If results are greater than 130 mmol/mol, the instrument will display results as > 130 mmol/mol. If your result is unusually high or low from your A1c goal set by the healthcare provider, immediately consult your healthcare professional.

POSSIBLE CAUSES OF ERROR AND INTERFERENCE

For information regarding troubleshooting and error messages, refer to the operator's guide for assistance.

NOTES FOR HEALTHCARE PROFESSIONALS

CALIBRATION

Instrument The DCA System is calibrated during manufacturing. Thereafter, the instrument performs self-checks during each assay. If the system detects an error, an error message displays.

Reagent Cartridges Reagent cartridges undergo thorough analysis and characterization during manufacturing. The DCA HbA1c test method is National Glycohemoglobin Standardization Program (NGSP) certified and is traceable to IFCC reference materials and test methods. The values for the calibration parameters are encoded onto the calibration card provided with each lot of reagent cartridges. This calibration information is transferred to the system by scanning the calibration card. The instrument can store calibrations for multiple lots. The system checks for a calibration curve for the lot in use when each reagent cartridge barcode containing the lot number and test name is scanned. If no calibration curve is found for the lot of reagent cartridges in use, the instrument prompts you to scan the calibration card.

PERFORMING QUALITY CONTROL

To ensure the quality of both the testing procedures and the patient results for hemoglobin A1c, the DCA System performs 48 optical, electronic, mechanical, and reagent systems checks during the course of each specimen assay. These checks include calibration verification during every test. If an assay or system error occurs during any individual measurement, the system automatically reports an error message, preventing the reporting of erroneous patient results.

Quality control (QC) testing confirms that your DCA System is working properly and providing reliable results. Only when controls are used routinely and the values are within acceptable ranges can accurate results be assured for patient samples.

Run control specimens under the following conditions:

- When using a new shipment of reagents
- When using a new lot number of reagents
- When training new operators
- Whenever test results are in doubt (for example, when results don't match the patient's clinical condition or symptoms)
- Every 30 days at sites with a low frequency of testing

Compare QC results to those listed as acceptable by the QC manufacturer. If QC results aren't acceptable, don't test patient samples until the problem is resolved. Repeat QC testing until the results are acceptable.

INTERPRETING RESULTS

If results are less than 4 mmol/mol, the instrument will display results as < 4 mmol/mol. These types of results are rare and may indicate that the sample contains substantial amounts of fetal hemoglobin; or that the patient may be suffering from hemolytic anemia or polycythemia. Record the result and consult a physician.

If results are greater than 130 mmol/mol, the instrument will display results as > 130 mmol/mol. All laboratory tests are subject to random error. If the test result is questionable or if clinical signs and symptoms appear inconsistent with test results, retest with a new sample or confirm the result using another method. Record the result and consult a physician.

LIMITATIONS AND KNOWN INTERFERENCES

The DCA HbA1c assay gives accurate and precise results over a range of total hemoglobin of 7–24 g/dL. Most patients will have hemoglobin concentrations within these values; however, patients with severe anemias may have hemoglobin concentrations < 7 g/dL, and patients with polycythemia may have hemoglobin concentrations > 24 g/dL.

Patients known to have these conditions should be assayed by a test employing a different assay principle if their hemoglobin concentrations are outside of the acceptable range.

Glycated hemoglobin F isn't measured by the DCA HbA1c assay. At levels of hemoglobin F < 10%, the DCA system accurately indicates the patient's glycemic control; however, at very high levels of hemoglobin F > 10%, the amount of HbA1c is lower than expected because a greater proportion of the glycated hemoglobin is in the form of glycated hemoglobin F. HbA1c results for such patients don't accurately indicate the patient's glycemic control and shouldn't be compared to published normal or abnormal values.

Conditions such as hemolytic anemia, polycythemia, homozygous HbS, and HbC, can result in decreased life span of the red blood cells, which causes HbA1c results to be lower than expected, regardless of the method used, and not to be related to glycemic control, when using published reference ranges.

Bilirubin up to a level of 30 mg/dL, doesn't interfere with this assay.

Triglycerides up to 1350 mg/dL in fresh whole blood, don't interfere with this assay.

Rheumatoid factor up to 1:33,000 titer, doesn't interfere with this assay.

Expected serum levels of the following drugs commonly prescribed to persons with diabetes don't interfere with this assay: **Diabinese, Orinase, Tolinase, Micronase, Dymelor, glipizide.**

EXPECTED VALUES

The expected normal range for HbA1c using the DCA HbA1c test was determined by assaying blood samples from 103 apparently healthy individuals (fasting blood glucose < 120 mg/dL). No significant differences in normal range were observed among males and females, geographical location, or age groups evaluated. The mean HbA1c value was (1 SD) 31 mmol/mol. The range was 22–48 mmol/mol. The 95% confidence limits (mean ± 2 SD) 23–39 mmol/mol. These values are similar to those reported in the literature.³

Depending on the assay methodology used, HbA1c is approximately 9–42 mmol/mol in non-diabetics, 42–64 mmol/mol in controlled diabetics and can be as much as 195 mmol/mol or higher in poorly controlled diabetics.⁴ However, each laboratory should determine normal ranges to conform to the population being tested.

PERFORMANCE CHARACTERISTICS

The precision and correlation data are results of studies conducted by the staff at separate physician offices. The statistical calculations were performed following Clinical Laboratory Standards Institute (CLSI) procedures.

Precision Results from studies conducted with multiple DCA HbA1c assays of 2 different commercially prepared blood controls were performed by 3 independent investigators. The imprecision of results using a commercially prepared normal or abnormal control was ≤ 1.2% with the normal control and ≤ 4.3% with the abnormal control (measured in %HbA1c results).

Correlation Clinical studies were conducted in which capillary blood results from the DCA 2000+ were compared with results obtained using an HPLC laboratory reference method. The slopes were found between 0.91 and 0.94 (when measured in %HbA1c).

SPECIFICITY

Effect of Hemoglobin Variants The antibody in the DCA HbA1c assay is specific for the first few amino acid residues of the glycated amino-terminus of the β-chain of hemoglobin A. Any glycated hemoglobin molecule having this same structure will be measured in the assay. Most glycated hemoglobin variants are immunoreactive in the DCA HbA1c assay (such as HbS1c, HbC1c, and HbE1c). The point mutations in these molecules occur at the 6 position of the β-chain (HbS and HbC) and at the 26 position of the β-chain (HbE). Thus, the point mutations in these variants don't affect the binding of the antibody used in the DCA HbA1c assay. The DCA reports HbA1c values that reflect the glycemic control of patients with these hemoglobinopathies.^{5,6}

Effect of Pre-HbA1c (Labile Fraction) The labile fraction (Schiff base attachment of glucose to HbA, or pre-HbA1c) doesn't affect the assay result because the antibody is specific for the stable ketoamine.⁵

Effect of Carbamylated Hemoglobin Carbamylated hemoglobin (elevated in patients with uremia) doesn't affect the assay result because the antibody is specific for the sugar moiety of HbA1c.^{7,8}

TECHNICAL ASSISTANCE

For customer support, please contact your local support provider or distributor.

REFERENCES

- Standard of Medical Care of Diabetes -DIABETES CARE, VOLUME 36, SUPPLEMENT 1, JANUARY 2013.
- Larsen, M. L., Horder, M., and Mogensen, E. F. Effect of long-term monitoring of glycosylated hemoglobin levels in insulin-dependent diabetes mellitus. *N. Eng. J. Med.* 323 (1990): 1021–1025.
- Nathan, D. M., Singer, D. E., Huxthal, K., and Goodson, J. D. The clinical information value of the glycosylated hemoglobin assay. *N. Eng. J. Med.* 310 (1984): 341–346.
- Goldstein, D. E., Little, R. R., and Wiedmeyer, H. M., et al. Glycated hemoglobin: Methodologies and clinical applications. *Clin. Chem.* 32 (1986): 864–870.
- Weykamp, C. W., Pender, T. J., Muskiet, F. J., and Van der Silk, W. Influence of hemoglobin variants and derivatives on glycohemoglobin determinations, as investigated by 102 laboratories using 16 methods. *Clin. Chem.* 39 (1993): 1717–1722.
- Eaton, S. E., Fielden, P., and Haisman, P. Glycated hemoglobin (HbA1c) measurements in subjects with hemoglobin variants, using the DCA 2000. *Ann. Clin. Biochem* 34 (1997): 205–207.
- Rose, A., Tongate, C., and Valdes, Jr., R. A hemoglobin A1c immunoassay method not affected by carbamylated hemoglobin. *Ann. Clin. Lab. Science* 25 (1995): 13–19.
- Weykamp, C. W., Miedema, K., de Haan, T., and Doelman, C. J. A. Carbamylated Hemoglobin Interference in Glycohemoglobin Assays. *Clin. Chem.* 45 (1999): 438–440.

Origin US

Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591-5097 USA

Siemens Healthcare Diagnostics
Manufacturing Ltd.
Chapel Lane
Swords, Co. Dublin, Ireland
www.siemens.com/poc

Siemens Healthineers
Headquarters
Siemens Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany
Phone: +49 91 31 84-0
siemens-healthineers.com

Date of Issuance 2019-07

© 2014–2019 Siemens Healthcare Diagnostics. All rights reserved.
DCA and DCA Vantage are trademarks of Siemens Healthcare Diagnostics.
All other trademarks are the property of their respective owners.

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Siemens Healthcare Diagnostics Inc. 511 Benedict Avenue Tarrytown, NY 10591-5097 USA