

INNOVANCE® D-Dimer Sample Diluent

INNOVANCE D-Dimer **DILUENT**



Revision bar indicates update to previous version.

Intended Use

INNOVANCE D-Dimer **DILUENT** is a buffer solution intended to be used in conjunction with the INNOVANCE® D-Dimer method, for the measurement of samples outside the initial measuring range.

Reagents

Reagent	Description	Storage	Stability
INNOVANCE® D-Dimer Sample Diluent INNOVANCE D-Dimer DILUENT	Ready to use liquid containing: • Buffers, preservatives	2–8 °C May be used up to the expiry date indicated on the label if stored unopened.	2–8 °C: once opened, 4 weeks ^a ; ≤ –18 °C: once opened, 4 weeks ^a

^a closed original vial

Warnings and Precautions

For *in-vitro* diagnostic use only.

For laboratory professional use.

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

Safety data sheets (MSDS/SDS) available on [siemens-healthineers.com/sds](https://www.siemens-healthineers.com/sds).



Danger! INNOVANCE D-Dimer **DILUENT**

Hazardous ingredient: Imidazole (0.332 % [w/w]).

H360D: May damage the unborn child.

P201: Obtain special instructions before use. **P280:** Wear protective gloves/protective clothing/eye protection/face protection. **P308 + P313:** IF exposed or concerned: Get medical advice/attention.

Contains sodium azide as a preservative. Sodium azide can react with copper or lead plumbing to form explosive metal azides. On disposal, flush reagents with a large volume of water to prevent the buildup of azides. Disposal into drain systems must be in compliance with prevailing regulatory requirements.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with all government requirements.

Preparing Reagents

- INNOVANCE D-Dimer **DILUENT** is ready for use. Prior to placing on the system: Mix carefully and avoid foam formation.
- INNOVANCE D-Dimer **DILUENT** is placed on the system prior to testing. Please refer to the Instructions for Use of INNOVANCE® D-Dimer for further details.

Procedure

Materials Provided

REF	Contents		
OPBR03	INNOVANCE® D-Dimer Sample Diluent INNOVANCE D-Dimer DILUENT	10 ×	5 mL

Materials Required but not Provided

Item	Description
REF OPBP03, OPBP07	INNOVANCE® D-Dimer and associated materials according to the INNOVANCE® D-Dimer Instructions for Use.
Coagulation analyzers ^b , such as:	<ul style="list-style-type: none"> • Atellica® COAG 360 System • BCS® XP System • SYSMEX CA-500/CA-600 series • SYSMEX CA-1500 System • SYSMEX CS-2000i/CS-2100i System • SYSMEX CS-2500 System • SYSMEX CS-5100 System

^b Availability of analyzers may vary by country.

Please note that the applications on other analyzers can be validated by the instrument manufacturer in accordance with the requirements of the REGULATION (EU) 2017/746 under their responsibility as long as the intended purpose and performance are not modified.

























- The sample dilution with INNOVANCE D-Dimer **DILUENT** is conducted automatically by the respective system/analyzer for samples with D-dimer concentrations within the extended measuring range.
- Please refer to the INNOVANCE® D-Dimer Instructions for Use and the respective Application Sheets.

Technical Assistance

For customer support, contact your local technical support provider or distributor.
siemens-healthineers.com

Definition of Symbols

The following symbols may appear on the product labeling:

	Do not reuse		Use By
	Batch Code		Catalogue Number
	Caution		Manufacturer
	Authorized representative in the European Community		Contains sufficient for <n> tests
	Biological Risks		<i>In Vitro</i> Diagnostic Medical Device
	Temperature Limitation		Consult instruction for Use
	Non-sterile		CE marking of conformity
	CE marking of conformity with notified body ID number. Notified body ID number can vary.		Contents
	Reconstitution volume		Level
	Keep away from sunlight and heat		Warning
	Danger		Prescription device (US only)
	Device Identification (UDI) barcode		REACH Authorization Number

Legal Information

Atellica, BCS and INNOVANCE are trademarks of Siemens Healthineers.

SYSMEX is a trademark of SYSMEX CORPORATION.

All other trademarks are the property of their respective owners.

© Siemens Healthineers, 2009–2021. All rights reserved.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestraße 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com



Siemens Healthcare Diagnostics Products GmbH

Emil-von-Behring-Str. 76
35041 Marburg
Germany
siemens-healthineers.com