

# Pathromtin® SL

C€0197

| Revision bar indicates update to previous version.

## Intended Use

Pathromtin® SL is an in vitro diagnostic reagent for the quantitative determination of activated partial thromboplastin time (APTT) as an aid to diagnosis, screening for hemostasis disorders and monitoring of unfractionated heparin in human sodium citrated plasma by means of automated, semi-automated and/or manual coagulometric methods.

For APTT testing no international reference preparation or method is available.

## Summary and Explanation

Pathromtin® SL is a liquid reagent based on plant phospholipids with silicon dioxide particles for plasma activation.

The APTT, a global screening procedure<sup>1,2</sup> used primarily to evaluate coagulation abnormalities in the intrinsic pathway, will also detect severe functional deficiencies in factors: FII, FV, FX, or fibrinogen. The APTT has also been widely advocated<sup>3-6</sup> as a means to monitor the effectiveness of unfractionated heparin therapy where the clotting time is prolonged in proportion to the level of heparin. In patients receiving oral anticoagulants, the circulating levels of factors: FII, FVII, FIX, and FX are depressed therefore the APTT can be expected to be prolonged. The presence of non-specific inhibitors, such as the lupus-like anticoagulant<sup>1,7</sup>, may prolong the APTT but this effect is variable and generally recognized as being related more to the nature of the APTT reagent employed.

In summary, the APTT is a clinically important screening test with wide applicability for the diagnosis of coagulant disorders and therapeutic monitoring of both, hemorrhagic and thrombotic disease<sup>8-10</sup>, and is specifically used

- In pre-surgery bleeding risk assessment
- In screening for bleeding disorders, e.g. in case of suspected intrinsic factor deficiency or inhibitors to intrinsic coagulation factors
- As an aid to diagnosis of lupus anticoagulants (inhibitors) in patients with thrombophilia
- For monitoring of therapy in patients receiving unfractionated heparin

Furthermore, Pathromtin® SL can be used in combination with the respective factor: FVIII, FIX, FXI or FXII deficient plasma for the quantification of the coagulation factors: FVIII, FIX, FXI and FXII.

## Principles of the Procedure

Incubation of plasma with the optimal quantity of phospholipids and a surface activator leads to activation of factors of the intrinsic coagulation system. The addition of calcium ions triggers the coagulation process; the time to formation of a fibrin clot is measured.

## Reagents

**Note:** Pathromtin® SL can be used on automated coagulation analyzers. Siemens Healthineers provides Reference Guides (Application Sheets) for several coagulation analyzers. The Reference Guides (Application Sheets) contain analyzer/assay specific handling and performance information

which may differ from that provided in these Instructions for Use. In this case, the information contained in the Reference Guides (Application Sheets) supersedes the information in these Instructions for Use. Please also consult the instruction manual of the instrument manufacturer!

Reagent	Description	Storage	Stability
Pathromtin® SL	Ready to use liquid containing: <ul style="list-style-type: none"> <li>• Silicon dioxide (1.2 g/L)</li> <li>• soy lecithin (0.25 g/L)</li> <li>• Sodium chloride</li> <li>• HEPES, pH 7.6</li> <li>• Preservative: <ul style="list-style-type: none"> <li>• Sodium azide (&lt; 1 g/L)</li> </ul> </li> </ul>	2–8 °C May be used up to the expiry date indicated on the label if stored unopened.	2–25 °C: once opened, 2 weeks <sup>a</sup>

<sup>a</sup> closed original vial

**Signs of expiry:** Deviations from the normal laboratory value in the determination of normal plasma or controls.

#### On-board stability

Information regarding on-board stability is specified in the Reference Guides (Application Sheets) for the different coagulation analyzers.

## Warnings and Precautions

For *in-vitro* diagnostic use only.

For laboratory professional use.

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

Safety data sheets (MSDS/SDS) available on [siemens-healthineers.com/sds](https://www.siemens-healthineers.com/sds).

Contains sodium azide as a preservative. Sodium azide can react with copper or lead plumbing to form explosive metal azides. On disposal, flush reagents with a large volume of water to prevent the buildup of azides. Disposal into drain systems must be in compliance with prevailing regulatory requirements.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with all government requirements.

Summary of Safety and Performance (SSP) is available in the European database on medical devices (see Eudamed public website: <https://ec.europa.eu/tools/eudamed>). In case Eudamed is not available, SSP can be delivered by Siemens Healthineers on request.

## Preparing Reagents

Pathromtin® SL must be gently inverted (5 to 8 times) to mix before first use.

Pathromtin® SL must be used at room temperature (15 to 25 °C).

Every 24 hours of use, the reagent must be inverted to resuspend any sediment.

CaCl<sub>2</sub> **SOLUTION** 0.025 mol/L: warm to 37 °C.

## Specimen Collection and Handling

### Collecting the Specimen

To obtain the plasma, carefully mix one part of 0.11 or 0.13 mol/L (3.2 % or 3.8 %) sodium citrate solution with nine parts of freshly collected patient blood, avoiding the formation of foam. It is recommended that blood specimens for plasma-based coagulation testing should be collected by venipuncture using a blood collection system that collects the specimen directly into a glass or plastic evacuated tube containing the appropriate additive. Immediately centrifuge for no less than 15 minutes at 1 500 × g, remove the supernatant plasma and keep at 15 to 25 °C until use in the test (max. 4 hours).

## Procedure

### Materials Provided

REF	Contents		
OQGS29	Pathromtin® SL	10 ×	5 mL
OQGS35	Pathromtin® SL	20 ×	5 mL

### Materials Required but not Provided

Item	Description
[REF] ORHO37	CaCl <sub>2</sub> [SOLUTION], Calcium Chloride Solution, (0.025 mol/L)
[REF] ORKE41	[CONTROL N], Control Plasma N, or
[REF] 291070	Dade® Ci-Trol® 1, or
[REF] B4244-10	Ci-Trol [CONTROL 1], Dade® Ci-Trol® Coagulation Control Level 1, as control for the normal range
[REF] OUPZ17	[CONTROL P], Control Plasma P (only for Atellica® COAG 360 System, BFT // Analyzer, BCS® XP System), or
[REF] 291071	Dade® Ci-Trol® 2, or
[REF] B4244-20	Ci-Trol [CONTROL 2], Dade® Ci-Trol® Coagulation Control Level 2, as control for the pathological/therapeutical range
[REF] 291072	Dade® Ci-Trol® 3, or
[REF] B4244-30	Ci-Trol [CONTROL 3], Dade® Ci-Trol® Coagulation Control Level 3, as control for the pathological/therapeutical range
[REF] B4224-50	Ci-Trol [HEPARIN] [CONTROL 1], Dade® Ci-Trol® Heparin Control, Low
[REF] B4224-60	Ci-Trol [HEPARIN] [CONTROL 2], Dade® Ci-Trol® Heparin Control, High
–	Pipettes for precise measurement of 0.1 mL
–	Plastic test tubes
Coagulation analyzers <sup>b</sup> , such as:	<ul style="list-style-type: none"> <li>• Atellica® COAG 360 System</li> <li>• BCS® XP System</li> <li>• BFT // Analyzer</li> <li>• SYSMEX CA-500/CA-600 series</li> <li>• SYSMEX CA-1500 System</li> <li>• SYSMEX CS-2000i/CS-2100i System</li> <li>• SYSMEX CS-2500 System</li> <li>• SYSMEX CS-5100 System</li> </ul>

<sup>b</sup> Availability of analyzers may vary by country.

Please note that the applications on other analyzers can be validated by the instrument manufacturer in accordance with the requirements of the REGULATION (EU) 2017/746 under their responsibility as long as the intended purpose and performance are not modified.

### Manual Testing

Pre-warm CaCl <sub>2</sub> [SOLUTION] at 37 °C	
Pipette into a test tube pre-warmed to 37 °C:	
Citrated plasma	100 µL

Pathromtin® SL	100 µL Mix well. Incubate at 37 °C for 2 minutes.
CaCl <sub>2</sub> <b>SOLUTION</b> (pre-warmed)	100 µL On addition of CaCl <sub>2</sub> <b>SOLUTION</b> start stop-watch or timer on the coagulation analyzer and determine the coagulation time.

## Monitoring of Unfractionated Heparin Therapy with APTT

When using the APTT for this purpose, the factors influencing the test should be kept in mind. General considerations are listed below.

- A. Time of collection is important since the *in-vivo* half-life of unfractionated heparin is approximately 1.5 hours<sup>5</sup>. When it is administered, it has an immediate anticoagulant effect but the degree of this effect decreases rapidly with time. This is especially apparent with intermittent single intravenous injections.
- B. The anticoagulant used for sample collection can alter test results.
- C. Platelet factor 4, a heparin neutralizing factor in platelet alpha-granules, can be released by platelet aggregation or damage. To prevent this occurrence *in-vitro*, the specimen should be collected with a minimum of trauma. Cold temperatures are known to induce platelet aggregation and release platelet factor 4; therefore, centrifugation at room temperature is recommended for heparin testing.
- D. Using APTT to monitor unfractionated heparin therapy is time-dependent. Delay in testing samples will result in prolonged APTT determinations. Therefore, it is imperative that the testing on all samples be performed as soon as possible.
- E. Increased contact activation times may result in prolonged APTT in plasma containing heparin. It is imperative that the optimal heating-activation time of the Pathromtin® SL-plasma mixture be rigidly standardized<sup>12</sup>.
- F. Different test systems (i.e. manual, photo-optical, etc.) will exhibit variable heparin sensitivity. Interchanging of test systems should be avoided.
- G. Baseline data on the APTT of each patient before the start of therapy should be established where feasible to determine the respective patient APTT as it relates to the reference range established for the test in that laboratory.
- H. Studies have shown variability in original estimates of the quality of unfractionated heparin from different sources and different manufacturers. *In-vivo* reactivity varies with the type of heparin administered, the metabolism of the individual and other coadministered medications<sup>5,12</sup>.
- I. Because the APTT can vary with technique, method, equipment, reagent lot and heparin used, each laboratory must establish its own therapeutic ranges, or verify them whenever one or more of the aforementioned variables is changed. This can be done by simultaneously determining the APTT and the heparin concentration for samples from patients receiving heparin therapy. A dose response curve can be calculated from the data using regression analysis, and the APTT range corresponding to a heparin concentration of 0.3 to 0.7 U/mL (by a factor Xa inhibition assay) can be derived<sup>5,13,4</sup>.

## Internal Quality Control

- Normal range: Dade® Ci-Trol® 1, Ci-Trol **CONTROL 1**, or **CONTROL N**
- Pathological range: Dade® Ci-Trol® 2, Ci-Trol **CONTROL 2**, or Dade® Ci-Trol® 3, Ci-Trol **CONTROL 3**, or **CONTROL P**
- Heparin monitoring: Ci-Trol **HEPARIN CONTROL 1**  
Ci-Trol **HEPARIN CONTROL 2**

Two controls (one in the normal range and one in the pathological/therapeutic range) must be measured at the start of the test run, after each change of reagent vial, and at least once during an 8-hour shift. The control material should be prepared and processed in the same manner as the patient samples. Each laboratory should establish its own confidence intervals for the controls. This interval is generally  $\pm 2$  to  $\pm 2.5$  standard deviations from the mean control value. If the control values are outside of the confidence interval, the controls, reagents and instrument must be checked. Before reporting

the patient values, it is recommended that all steps should be documented that were taken to identify and rectify the problem. New control ranges should be defined for each new lot of reagents or controls.

## Results

The result is given in seconds.

## Limitations

It should be noted that APTT testing may be affected by several commonly administered drugs. Therapeutic doses of hirudin or other direct thrombin inhibitors may prolong clotting times<sup>18,19</sup>. Lipoglycopeptide antibacterial drugs (such as oritavancin or telavancin) may interfere with APTT based assays. Consult Instructions for Use of respective drugs<sup>20,21</sup>.

Decrease in time of APTT determination in conjugated estrogen therapy in males and oral contraceptive administration in females has been reported<sup>22,23</sup>. Increase in the APTT has been seen in diphenylhydantoin, heparin, warfarin, naloxone and radiographic agent administration<sup>24</sup>. Heparin or low molecular weight Heparin (LMWH) is used for the treatment of thromboembolic diseases followed by prophylaxis of recurrence using oral anticoagulants. The effect of Heparin is monitored by the activated partial thromboplastin time (APTT), the thrombin time (TT) or a dedicated Heparin assay. The discontinuation of heparin depends on whether the therapeutic range of the oral anticoagulation is achieved. Heparinase (Enzyme Commission number EC 4.2.2.7) has been shown to eliminate heparin from plasma<sup>25</sup>.

The APTT assay is a functional test which screens for global coagulation disorders of the endogenous coagulation system. It is common knowledge that a low concentration of coagulation factors: FII, FV, FX and high concentration of fibrin(ogen) degradation products also have influence on the assay<sup>26</sup>. The condition of the specimen may affect results. Hemolyzed, icteric and /or lipemic samples might either interfere chemically, biologically or with the physical properties of the detection principle of instruments<sup>27</sup>.

The latter is particularly true for optical instrumentation measurements of the APTT. Preactivation of the sample due to poor phlebotomy technique can lead to false results. Referencing the Sample Handling and Collection section of CLSI Document H21-A5<sup>11</sup> it is recommended that test tubes with non-wettable surfaces (e.g., plastic) be used rather than glass test tubes.

Siemens Healthineers has validated use of these reagents on various analyzers to optimize product performance and meet product specifications. Please note that the applications on other analyzers can be validated by the instrument manufacturer in accordance with the requirements of the REGULATION (EU) 2017/746 under their responsibility as long as the intended purpose and performance are not modified. User defined modifications are not supported by Siemens Healthineers as they may affect performance of the system and assay results. It is the responsibility of the user to validate modifications to these instructions or use of the reagents on analyzers other than those included in Siemens Healthineers Application Sheets or these Instructions for Use.

Results of this test should always be interpreted in conjunction with the patient's medical history, clinical presentation and other findings.

## Expected Values

For the purpose of APTT screening, results in seconds should be related to the reference interval for APTT testing in each laboratory.

Reference intervals vary from laboratory to laboratory depending on the population served and the technique, method, equipment and reagent lot used. Therefore, each laboratory must establish its own reference intervals or verify them whenever one or more of the aforementioned variables are changed.

In a study of ostensibly healthy individuals using a specific lot of Pathromtin® SL, the following values were obtained:

	n	Median [s]	90 % Reference Interval	
			5 <sup>th</sup> Percentile [s]	95 <sup>th</sup> Percentile [s]
SYSMEX CA-1500 System	100	30.9	26.4	37.5
BCS® System	111	30.2	25.9	36.6

Reference ranges for other populations such as pediatric groups should also be established where warranted.

**Note:** CLSI Document C28-A3c (cited in H47-A)<sup>29</sup>, states that a parametric approach (mean  $\pm$  2 SD) can be applied. The assumption of this approach (Gaussian normal distribution) must however be checked.

## Performance Characteristics

### Measuring Range

The measuring range depends on the individual application of the assay due to instrument related conditions. Application specific performance data are listed in the respective Reference Guides of the instruments.

### Sensitivity

For characterization of an APTT reagent, determination of heparin sensitivity, factor sensitivity, and Lupus sensitivity is important. These are reagent characteristics which are independent from the analyser used.

#### Factor Sensitivity

According to CLSI H47-A2, the APTT reagent/instrument combination used should provide abnormally prolonged results for plasmas that have less than 30 % factor activity of the coagulation factors: FVIII, FIX and FXI. CLSI H47-A2 recommends to determine sensitivity levels by serial dilution of normal plasma into deficient plasma. Sensitivity levels determined by this method should ideally be within 30 and 45 %. However, the factor sensitivity levels determined by this method is strongly dependent on deficient plasma used<sup>28</sup>.

#### Heparin Sensitivity

For determination of heparin sensitivity 167 samples from patients receiving unfractionated heparin were analyzed. By comparing the APTT ratio obtained with Pathromtin® SL and the Heparin activity measured with Berichrom Heparin a correlation coefficient  $r = 0.709$  was observed. Each individual laboratory hospital should determine its own therapeutic heparin range using the ex vivo method according the CLSI guidelines (H57-A2)<sup>29</sup>.

#### Lupus Sensitivity

By testing a panel of LA-negative donors from ostensibly healthy donors and 97 commercially available frozen plasma samples from LA-positive patients (as determined by dRVVT, PT and thrombin time) the following sensitivity was observed:

##### Sensitivity of Pathromtin® SL<sup>30</sup>

APTT Reagent	Sensitivity (%) 95 <sup>th</sup> cut-off	Sensitivity (%) 99 <sup>th</sup> cut-off
Pathromtin® SL	72.2	59.8

### Precision

Using opto-densitometric detection, precision studies were run over a 5 day period, twice per day ( $n = 8$  replicates per day), to total  $n = 40$  for normal and pathological control plasmas as well as a heparin plasma pool. Intra-assay precision was calculated from the individual  $n = 4$  precision values for the daily runs over the 5 days. Intra-assay precision ranged from 0.6 to 2.0 % CV, while the inter-assay precision ranged from 0.3 to 2.8 % CV.

Other system specific results are given in the respective Reference Guides (Application Sheets).

The reproducibility was assessed by Siemens Healthineers for APTT with Pathromtin® SL based on publicly available proficiency testing information in 2018/2019. The overall reproducibility median CV% was found to be <5 % including lot, instrument, laboratory and operator variability factors.

## Method Comparison

Results of a comparison of APTT determinations using Pathromtin® SL Reagent and another commercially available APTT reagent gave a correlation coefficient of 0.96 and a y-intercept of 2.1 s, and a slope of 0.99.

## Technical Assistance

For customer support, contact your local technical support provider or distributor.  
siemens-healthineers.com

## Current Version of Application Sheets

Pathromtin® SL can be used in combination with various automated coagulation analyzers. Siemens Healthineers provides Reference Guides/Application Sheets for the coagulation analyzers listed in section "Materials Required but not Provided", page 3 under the dedicated link below:  
siemens-healthineers.com/rg

As Siemens Healthineers continuously monitors the product performance and safety, the users are required to ensure that they work with the correct revision of the instructions for the product lots in use. Please periodically review the availability of new electronic labeling revisions to ensure safe use of the product.

The IFU version number is visible on each product box label. Siemens Healthineers ensures that all products lots bearing the same IFU version number are compatible with the electronic labeling provided via siemens-healthineers.com/eIFU.

## References

























1. Kamal AH, Tefferi A, Pruthi RK. How to interpret and pursue an abnormal prothrombin time, activated partial thromboplastin time, and bleeding time in adults. *Mayo Clin Proc.* 2007;82:864-73.
2. Leung LLK. Perioperative evaluation of bleeding diathesis. *Hematology Am Soc Hematol Educ Program.* 2006;457-61.
3. Eikelboom JW, Hirsh J. Monitoring unfractionated heparin with the aPTT: Time for a fresh look. *Thromb Haemost.* 2006;96:547-52.
4. Bates SM, Johnston M, Hirsh J, et al. Use of a fixed activated partial thromboplastin time ratio to establish a therapeutic range for unfractionated heparin. *Arch Intern Med.* 2001;161:385-91.
5. Hirsh J, Raschke R. Heparin and low-molecular-weight heparin. The seventh ACCP Conference on Antithrombotic and Thrombolytic Therapy. *Chest.* 2004;126:188S-203S.
6. Marlar RA, Clement B, Gausman J. Activated partial thromboplastin time monitoring of unfractionated heparin therapy: Issues and recommendations. *Semin Thromb Hemost.* 2017;43:253-60.
7. Teruya J, West AG, Suell MN. Lupus Anticoagulant Assays. Questions answered and to be answered. *Arch Pathol Lab Med.* 2007;131:885-9.
8. Fritsma GA, Dembitzer FR, Randhawa A, et al. Recommendations for appropriate activated partial thromboplastin time reagent selection and utilization. *Am J Clin Pathol.* 2012;137:904-8.
9. Kershaw G. Performance of Activated Partial Thromboplastin Time (APTT): Determining Reagent Sensitivity to Factor Deficiencies, Heparin, and Lupus Anticoagulants. *Methods Mol Biol.* 2017;1646:75-83.
10. Levy JH, Szlam F, Wolberg AS, et al. Clinical use of the activated partial thromboplastin time and prothrombin time for screening: a review of the literature and current guidelines for testing. *Clin Lab Med.* 2014;34:453-77.
11. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline – Fifth Edition. CLSI document **H21-A5**. Wayne, PA: Clinical and Laboratory Standards Institute; 2008.
12. Hattersley PG, Hayse D. The effect of increased contact activation time on the activated partial thromboplastin time. *Am J Clin Pathol.* 1976;66:479-82.

13. Olson JD, Arkin CF, Brandt JT, et al. College of American Pathologists Conference XXXI on laboratory monitoring of anticoagulant therapy: laboratory monitoring of unfractionated heparin therapy. *Arch Pathol Lab Med*. 1998;122:782-98.
14. Poller L, Thomson JM. The activated partial thromboplastin time (aPTT) In: Jespersen J, Bertina RR, Haverkate F, editors. *A Manual of Laboratory Techniques*. Hingham, MA: Kluwer Academic Publishers; 1992. p. 35-40.
15. Brack MJ, More RS, Hubner PJ, et al. The effect of low dose nitroglycerin on plasma heparin concentrations and activated partial thromboplastin times. *Blood Coagul Fibrinolysis*. 1993;4:183-6.
16. Harder S, Graff J, Klinkhardt U, et al. Transition from argatroban to oral anticoagulation with phenprocoumon or acenocoumarol: effects on prothrombin time, activated partial thromboplastin time, and Ecarin Clotting Time. *Thromb Haemost*. 2004;91:1137-45.
17. Boogen C, Niederau C, Reinauer H. Assessment of the influence of r-hirudin on coagulometric factor assays [Abstract]. *Ann Hematol*. 1999;78 Suppl:73.
18. Greinacher A, Warkentin TE. The direct thrombin inhibitor hirudin. *Thromb Haemost*. 2008;99:819-29.
19. Young DS, editor. *Effects of Drugs on Clinical Laboratory Tests*, 5th ed. Washington: AACC Press; 2000.
20. ORBACTIV (The Medicines Co.), Manufacturer information. Labels for NDA 206334 [cited 2016 Oct 17]. Available from: <https://www.accessdata.fda.gov/scripts/cder/daf/index.cfm>.
21. VIBATIV (Theravance Inc), Manufacturer information. Labels for NDA 022407 [cited 2013 Jun 21]. Available from: <https://www.accessdata.fda.gov/scripts/cder/daf/index.cfm>.
22. Prasad RN, Koh SC, Viegas OA, et al. Effects on hemostasis after two-year use of low dose combined oral contraceptives with gestodene or levonorgestrel. *Clin Appl Thromb Hemost*. 1999;5:60-70.
23. Ambrus JL, Schimert G, Lajos TZ, et al. Effect of antifibrinolytic agents and estrogens on blood loss and blood coagulation factors during open heart surgery. *J Med*. 1971;2:65-81.
24. Solomon GE, Hilgartner MW, Kutt H. Coagulation defects caused by diphenylhydantoin. *Neurology*. 1972;22:1165-71.
25. Harenberg J, Reichel T, Malsch R, et al. Multicentric evaluation of heparinase on aPTT, thrombin clotting time and a new PT reagent based on recombinant human tissue factor. *Blood Coagul Fibrinolysis*. 1996;7:453-458.
26. Thomas L, Chapter 16 in: *Clinical Laboratory Diagnostics (Labor und Diagnose)*. 8th ed. TH Books; 2016.
27. Lippi G, Franchini M, Montagnana M, et al. Quality and reliability of routine coagulation testing: can we trust that sample? *Blood Coagul Fibrinolysis*. 2006;17:513-9.
28. Lawrie AS, Kitchen S, Efthymiou M, et al. Determination of APTT factor sensitivity- the misleading guideline. *Int J Lab Hematol*. 2013;35:652-657.
29. NCCLS. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition. NCCLS document **C28-A2** [ISBN 1-56238-406-6]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA 2000
30. Wissel T, Wagner C. Lupus anticoagulants sensitivity of Siemens APTT reagents. *J Thromb Haemost*. 2013;Suppl.2:1167.



## Definition of Symbols

The following symbols may appear on the product labeling:

	Do not reuse		Use By
	Batch Code		Catalogue Number
	Caution		Manufacturer
	Authorized representative in the European Community		Contains sufficient for <n> tests
	Biological Risks		<i>In Vitro</i> Diagnostic Medical Device
	Temperature Limitation		Consult instruction for Use
	Non-sterile		CE marking of conformity
	CE marking of conformity with notified body ID number. Notified body ID number can vary.		Contents
	Reconstitution volume		Level
	Keep away from sunlight and heat		Warning
	Danger		Prescription device (US only)
	Device Identification (UDI) barcode		REACH Authorization Number

## Legal Information

Atellica, BCS, Ci-Trol, Dade and Pathromtin are trademarks of Siemens Healthineers.

SYSTEMEX is a trademark of SYSTEMEX CORPORATION.

All other trademarks are the property of their respective owners.

© Siemens Healthineers, 2010–2021. All rights reserved.

### Siemens Healthineers Headquarters

Siemens Healthcare GmbH  
Henkestraße 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
siemens-healthineers.com



### Siemens Healthcare Diagnostics Products GmbH

Emil-von-Behring-Str. 76  
35041 Marburg  
Germany  
siemens-healthineers.com