IMMULITE® Intact PTH Control Module

English

The IMMULITE® Intact PTH Control Module is an assayed, bi-level intact PTH control intended for use with the IMMULITE/IMMULITE 1000 (LKPP), IMMULITE 2000 (L2KPP and IMMULITE/IMMULITE 1000 Turbo (LSKPT) Intact PTH assays. It is intended as an aid in monitoring day-to-day assay performance.

Module Catalog Number: **LPHCM**

Control Lot: **0067**

Control Expiration: **2024-02**

Materials Supplied

Control 1 and 2 (PHCO1, PHCO2)

Six vials (three vials each) containing lyophilized synthetic human intact PTH in a buffered matrix, with preservative.

At most 30 minutes before use: reconstitute each vial with **2.0 mL** distilled or deionized water. Mix by *gentle* swirling or inversion until the lyophilized material is fully dissolved. After reconstitution, *place on ice immediately and use as soon as possible. **Discard after use.** (Do not aliquot and freeze.)*

The controls are intended to be assayed as unknowns, in the same manner as patient samples, in the context of an internal quality control program. Each control should be identified by its lot number (read from the vial) and the date it was opened/reconstituted.

Warnings and Precautions

For *in vitro* diagnostic use.

Do **not** use the control if it does not reconstitute completely, or if the reconstituted solution is not clear.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Sodium azide, at concentrations less than 0.1 g/dL, has been added as a preservative.

On disposal, flush with large volumes of water to prevent the buildup of potentially explosive metal azides in lead and copper plumbing.

Procedure

Assay in the same manner as patient samples, according to the procedures provided in the appropriate package inserts.

Target Values, pg/mL

IMMULITE/IMMULITE 1000 Intact PTH (LKPP)			
Level	Mean	SD	2SD Range
PHCO10067	58	5.8	46–70
PHCO20067	318	25	268–368
IMMULITE 2000 Intact PTH (L2KPP)			
Level	Mean	SD	2SD Range
PHCO10067	61	4.5	52–70
PHCO20067	340	32	276–404
IMMULITE/IMMULITE 1000 Turbo Intact PTH (LSKPT)			
Level	Mean	SD	2SD Range
PHCO10067	59	6.4	46–72
PHCO20067	320	31	258–382

Limitations

The control values should be used as a guide in evaluating performance. The control intervals and ranges should be adapted to each laboratory's individual requirements. Values obtained should fall within the established range. Each laboratory should establish corrective measures if individual values fall outside the range. Follow the applicable government regulations and local guidelines for quality control.

Technical Assistance

For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

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IMMULITE is a trademark of Siemens Healthcare Diagnostics.
© 2022 Siemens Healthcare Diagnostics. All rights reserved.