SIEMENS
Healthineers

IVD

CE

EN

INTERLOC USE

Chek-Stix UA Diptube Liquid QC is a ready-to-use liquid product intended for use in Near Patient Point-of-Care and Central Laboratory settings.

It is a control for qualitative and semi-quantitative procedures used in the biochemical and chemical determinations of human urine, for use only on the CLINITEK Status+ and CLINITEK Advanta® instruments.

SUMMARY AND EXPLANATION

Chek-Stix UA Diptube Liquid QC contains measured amounts of chemicals, stabilized human red cells and organic particles to simulate reagent strips used as a control for physical and chemical tests routinely performed in urinalysis. The Level 2 Positive Liquid QC contains <0.1% sodium citrate.

Kit Contents include: Level 1 Negative Urine Control and Level 2 Positive Urine Control.

STABILITY & STORAGE

Chek-Stix UA Diptube Liquid QC should be kept at 2-8°C when not in use. When stored at 2-8°C, the Control is stable for until the expiration date on the product labeling. The control can be used for Chemical analysis when stored at room temperature (20-25°C) for up to 30 days. After opening, store the control at 2-8°C until the expiration date stated on the label when stored at room temperature (2-8°C) between uses. Label the bottle with the date it was originally brought to room temperature to ensure proper performance. **DO NOT FREEZE.**

RISK AND SAFETY

WARNING!

H317 - May cause an allergic skin reaction.

P280 - Wash hands thoroughly after handling. Wear protective clothing, protective gloves, eye protection/face protection.

P333+P313 - If on skin or eyes: Wash with water.

P333+P313+P335 - If it irritates or causes rash: Get medical advice/attention.

P362+P364 - Take off immediately and wash with water.

Contains 2-methyldisothiazolin-3(2H)-one. (Level 2 positive).

Safety Data Sheets (MSDS/SDS) are available on seimens-healthineers.com/poc.

PRECAUTIONS

Because no method can offer complete assurance that HIV, HBsAg, HCV or other infectious agents are absent, it is recommended that human serum-based products be handled with the same precautions as the patient specimens.

Always follow your facility guidelines when working with biosafety materials.

AVAILABILITY

Chek-Stix UA Diptube Liquid QC is available in two different levels, Level 1 Negative and Level 2 Positive providing the laboratory a means of controlling and verifying accuracy and range of clinically significant values.

MATERIAL PROVIDED

• Chek-Stix UA Diptube Liquid QC, Level 1 Negative, 3 x 15ml.

• Chek-Stix UA Diptube Liquid QC, Level 2 Positive, 3 x 15ml.

• Value Sheet

• Special Materials Required (Not Provided)

• Test Tube

DIRECTIONS FOR USE

Standardized Chemical Tests

1. Remove the Chek-Stix UA Diptube Liquid QC from the refrigerator. Gently swirl the control to re-suspend any sediment and dispense an aliquot of no less than 1 mL into a test tube. Promptly recap the diptube and return it to the refrigerator.

2. Allow the Chek-Stix UA Diptube Liquid QC to reach room temperature (20-25°C) prior to testing.

3. Follow the CLINITEK Status+ , CLINITEK Advanta Operator or Reagent strip instructions for steps on performing a QC test.

NOTE: Allow no more than 20 drops into a single tube for the entry of the test strip.

A/C Testing

1. Follow steps 1-2 above.

2. Refer to the CLINITEK Status+ Operators Guide steps for performing a CLINITEK® IGG Pregnancy QC Test.

REAGENTS

The reagents used for quality control testing.

Standardized Chemical Tests

1. Remove the Chek-Stix UA Diptube Liquid QC from the refrigerator. Gently swirl the control to re-suspend any sediment and dispense an aliquot of no less than 1 mL into a test tube. Promptly recap the diptube and return it to the refrigerator.

2. Allow the Chek-Stix UA Diptube Liquid QC to reach room temperature (20-25°C) prior to testing.

3. Follow the CLINITEK Status+ , CLINITEK Advanta Operator or Reagent strip instructions for steps on performing a QC test.

NOTICE: Allow no more than 20 drops into a single tube for the entry of the test strip.

PROBLEMS/DIAGNOSTIC

If discrepancies arise from the expected ranges on the test specific insert, we recommend the following:

• Refer to user guides of the CLINITEK Status+ , CLINITEK Advanta, or your professional lot test reagent strip instructions.

• Immense the reagent in Chek-Stix UA Diptube Liquid QC and test a patient sample in accordance with the manufacturer's instructions.

• If the values remain beyond acceptable, try a different test, if possible, a different lot number of strips.

• If the discrepancy is in an instrument-generated value, check the instrument and check its calibration.

SYMBOLLEGEND

Consult Instructions for Use

IVD In Vitro Diagnostic Medical Device

Temperature Limits

LOT Lot Number

Caution, Consult Accompanying Documents

EC REP European Authorized Representative

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• Chek-Stix UA Diptube Liquid QC, Level 1 Negative, 3 x 15 ml.

• Chek-Stix UA Diptube Liquid QC, Level 2 Positive, 3 x 15 ml.

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