AlaTOP[®] Allergy Screen Control

English

The AlaTOP[®] Allergy Screen Controls are assayed, bi-level controls intended for use with the IMMULITE[®] 2000 AlaTOP Allergy Screen assay. They are intended as an aid in monitoring day-to-day assay performance.

Catalog Number: **L2ATC1–2**

Control Lot: **0135**

Control Expiration: **2023-07**

Materials Supplied

Alatop Allergy Screen Negative Control (L2ATC1)

2 mL IgE-free human serum. Stable at 2–8°C for 4 months after opening, or until the expiration date marked on the label.

AlaTOP Allergy Screen Positive Control (L2ATC2)

2 mL human serum with IgE specific for the IMMULITE 2000 AlaTOP Allergy Screen assay. Stable at 2–8°C for 4 months after opening, or until the expiration date marked on the label.

Do not freeze. Do not use if signs of microbiological contamination, such as a cloudy appearance, have occurred.

The controls are intended to be assayed as unknowns, in the same manner as patient samples, in the context of an internal quality control program. Each control should be identified by its lot number (read from the vial) and the date it was opened/reconstituted.

Warnings and Precautions

For *in vitro* diagnostic use.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Procedure

Assay in the same manner as patient samples, according to the procedures provided in the appropriate package inserts.

Target Values, Ratio

IMMULITE 2000 AlaTOP Allergy Screen (L2KAT)			
Level	Mean	SD	Ratio Range
L2ATC10135	0.45	0.225	0.0–0.90
L2ATC20135	3.3	0.33	2.6–4.0

Limitations

The control values should be used as a guide in evaluating performance. The control intervals and ranges should be adapted to each laboratory's individual requirements. Values obtained should fall within the established range. Each laboratory should establish corrective measures if individual values fall outside the range. Follow the applicable government regulations and local guidelines for quality control.

Technical Assistance

Available outside the United States only. For technical assistance, contact your National Distributor.

www.siemens.com/diagnostics

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

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IMMULITE and AlaTOP are trademarks of Siemens Healthcare Diagnostics.
All rights reserved.