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1. Austin L, Heath H. Calcitonin: physiology and pathophysiology. N Engl J Med. 1981; 304(5):269–278.
2. Zaidi M, Moonga BS, Bevis PJ, et al. Expression and function of the calcitonin gene products. Vitam Horm. 1991;46:87–164.
3. Faour O, Gilloteaux J. Calcitonin: Survey of new anatomy data to pathology and therapeutic aspects. Translational Research in

Anatomy. 2017;6:4–15.

4. Inzerillo AM, Ziadi M, Huang CL. Calcitonin: physiological actions and clinical applications. *J Pediatr Endocrinol Metab*. 2004; 17(7):931–940.
5. Boron WF, Boulpaep EL. *Medical Physiology: A Cellular And Molecular Approach*. Philadelphia, PA: Saunders; 2004.
6. Guesgen C, Willms A, Zwad A, et al. Investigation of factors potentially influencing calcitonin levels in the screening and follow-up for medullary thyroid carcinoma: a cautionary note. *BMC Clinical Pathology*. 2013; Nov 4; 13(1):27:1–8.
7. Society for Endocrinology. You and your hormones; calcitonin. 2011. <http://www.yourhormones.info/hormones/calcitonin.aspx>. Reviewed Jan 2015. Accessed June 20, 2017.
8. Hahm JR, Lee MS, Min YK, et al. Routine measurement of serum calcitonin is useful for early detection of medullary thyroid carcinoma in patients with nodular thyroid diseases. *Thyroid*. 2001; 11(1):73–80.
9. Pacini F, Fontanelli M, Fugazzola L, et al. Routine measurement of serum calcitonin in nodular thyroid diseases allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma. *J Clin Endocrinol Metab*. 1994 Apr;78(4):826–829.
10. Hillyard C, Oscier DG, Foa R, et al. Immunoreactive calcitonin in leukaemia. *Br Med J*. 1979 Dec 1;2(6202):1392–1393.
11. Hillyard C, Stevenson J. Calcitonin as a tumour marker. *Ligand Rev*. 1980(Spring); 2(1): 31–33.
12. Milhaud G. Ectopic secretion of calcitonin. In: *Calcitonin 1980; Chemistry, Physiology, Pharmacology, and Clinical Aspects*. Pecile A, ed. Amsterdam: Excerpta Medica. 1981: 154–169.
13. Mulder H, Hackeng WHL. Ectopic secretion of calcitonin. *Acta Med Scand*. 1978; 204: 253–256.
14. Schwartz KE, Wolfson AR, Forster B, et al. Calcitonin in nonthyroidal cancer. *J Clin Endocrinol Metab*. 1979 Sep;49(3):438–444.
15. Wallach S, Royston I, Taetle R, et al. Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung. *J Clin Endocrinol Metab*. 1981; Sep;53(3):602–606.
16. Mulder H, Silberbusch J, Hackeng WH, et al. Hypercalcitoninaemia in patients with chronic inflammatory disease. *Neth J Med*. 1980;23(3):129–131.
17. Toffaletti J. Calcitonin. *Clin Chem News*. 1984 (July);10(7):20–21.
18. Williams ED, Brown CL, Doniach I. Pathological and clinical findings in a series of 67 cases of medullary carcinoma of the thyroid. *J Clin Pathol*. 1966; Mar;19(2):103–113.
19. Hadoux J, Schlumberger M. Chemotherapy and tyrosine-kinase inhibitors for medullary thyroid cancer. *Best Pract Res Clin Endocrinol Metab*. 2017 Jun;31(3):335–347.
20. Leboulleux S, Baudin E, Travagli JP, et al. Medullary thyroid carcinoma. *Clin Endocrinol (Oxf)*. 2004 Sep;61(3):229–310.
21. Elisei R, Pinchera A. Advances in the follow-up of differentiated or medullary thyroid cancer. *Nat Rev Endocrinol*. 2012;8:466–475.
22. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections;

Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.

23. Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
24. Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
25. Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
26. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45(7):942–956.
27. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38(9):1737–1742.
28. Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
29. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.
30. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
31. Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
32. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

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Siemens Healthcare Diagnostics Inc.

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