

Squamous Cell Carcinoma (SCC)

| | | | |
|----------------------------------------------|------------------------------------------------------------|------------|-------------------------|
| Current Revision and Date^a | Rev. 02, 2022-06 | | |
| Product Name | ADVIA Centaur Squamous Cell Carcinoma (SCC) | REF | 11354594 (100 tests) |
| Abbreviated Product Name | ADVIA Centaur SCC | | |
| Test Name/ID | SCC | | |
| Systems | ADVIA Centaur CP system | | |
| Materials Required but Not Provided | ADVIA Centaur Squamous Cell Carcinoma Calibrator (SCC CAL) | REF | 11553960 |
| | ADVIA Centaur Wash 1 (2 × 1500 mL) | REF | 01137199 (112351) |
| | ADVIA Centaur Wash 1 (2 × 2500 mL) | REF | 03773025 |
| Optional Materials | ADVIA Centaur Tumor Marker Quality Control (TM QC) | REF | 11538186 |
| | ADVIA Centaur Multi-Diluent 13 | REF | 10492364 |
| Specimen Types | Serum, EDTA plasma | | |
| Sample Volume | 50 µL | | |
| Measuring Interval | 0.60–70.00 ng/mL (0.60–70.00 µg/L) | | |

^a A vertical bar in the page margin indicates technical content that differs from the previous version.



WARNING

The concentration of SCC in a given specimen, as determined by assays from different manufacturers, can vary due to differences in assay methods and reagent specificity. The results reported by the laboratory to the physician must include the identity of the SCC assay used. Values obtained with different assay methods cannot be used interchangeably. If, in the course of monitoring a patient, the assay method used for determining serial levels of SCC is changed, the laboratory must perform additional testing to confirm baseline values.

WARNING

SCC reactive determinants are shed naturally in saliva and other bodily fluids. Contamination of the samples or assay materials may cause falsely elevated SCC assay values.

Intended Use

The ADVIA Centaur® Squamous Cell Carcinoma (SCC) assay is for *in vitro* diagnostic use in the quantitative measurement of squamous cell carcinoma (SCC) antigen in human serum and plasma (EDTA) using the ADVIA Centaur® CP system.

The assay is used as an aid in the management of patients with squamous cell carcinoma. The measurement of SCC antigen, in conjunction with other clinical and laboratory findings, is used as an aid in monitoring disease progression during the course of disease and treatment in squamous cell carcinoma in lung, cervical, and head and neck cancer patients. Serial testing for patient SCC assay values should be used in conjunction with other clinical methods used for monitoring lung cancer.

Summary and Explanation

Squamous Cell Carcinoma Antigen (SCC Ag) is a 42 kDa subfraction of tumor-antigen-4 (TA-4), a glycoprotein with a molecular weight of approximately 48 kDa and multiple subfractions. The SCC Ag subfraction is a tumor marker extracted from squamous cell carcinomas of the uterine cervix.¹ The gene for SCC Ag encodes two antigens, SCC Ag1 and SCC Ag2, both detectable in serum.^{2,3} These antigens, present in the cytosol, are released into the bloodstream and may be elevated in patients with squamous cell carcinoma. SCC Ag1 and SCC Ag2 are co-expressed in the lung, head and neck, as well as the squamous epithelium of tongue, tonsils, esophagus, uterine cervix, and vagina.⁴

SCC Ag has been studied in squamous cell malignancies of the lung, esophagus, head and neck, anal canal, and skin. It has been reported that rising antigen in serial determinations may indicate disease recurrence. Measurement of SCC Ag levels has been used to monitor residual disease post treatment as well as response to therapy. More advanced cancer stages are associated with higher SCC Ag levels.⁵⁻¹¹

SCC Ag has specifically been reported as a biomarker for non-small cell lung cancer (NSCLC), primarily of the squamous cell carcinoma type. In NSCLC, SCC Ag has been used to monitor disease recurrence and residual disease following treatment and response to therapy.

SCC Ag is often used as a biomarker for cervical cancer and is recognized as the marker of choice for the follow-up of cervical cancer according to the European Group of Tumor Markers guidelines.¹² Elevated SCC Ag is associated with radiotherapy resistance and studies have demonstrated that the rate of SCC Ag reduction during radiotherapy can predict tumor response after treatment.¹³ Moreover, failure of SCC Ag levels to normalize posttreatment have been shown to predict tumor relapse with a high specificity.

Squamous cell-based cancers represent a significant number of head and neck cancers originating from the mucosal lining (epithelium) of these regions. Various studies have shown that SCC Ag is a significant independent predictor of disease-free survival and pretreatment levels are an independent prognostic indicator in these patients.

Principles of the Procedure

The ADVIA Centaur SCC assay is a fully automated, 2-site sandwich immunoassay using direct chemiluminescent technology. The Solid Phase contains magnetic microparticles coated with anti-SCC107 mouse monoclonal antibody. The Lite Reagent consists of acridinium ester-labeled anti-SCC140 mouse monoclonal antibody.

A direct relationship exists between the amount of SCC present in the patient sample and the amount of relative light units (RLUs) detected by the system.

Reagents

| Material Description | Storage | Stability |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------------------------|
| ADVIA Centaur SCC ReadyPack® primary reagent pack^{a,b} Lite Reagent 10.0 mL/reagent pack Anti-SCC monoclonal antibody covalently labeled with acridinium ester (~0.75 µg/mL) in buffer containing bovine serum albumin; surfactant; preservatives | Unopened at 2–8° C | Until expiration date on product |
| Solid Phase 17.5 mL/reagent pack Anti-SCC monoclonal antibody (< 0.001%) covalently coupled to magnetic microparticle in buffer containing bovine serum albumin; surfactant; preservatives | Onboard | 30 days |
| ADVIA Centaur Multi-Diluent 13 ReadyPack ancillary reagent pack^{a,b,c} 2 x 10.0 mL/pack Buffer; surfactant; sodium azide (< 0.1%) | Unopened at 2–8° C | Until expiration date on product |
| | Onboard | 28 days |
| ADVIA Centaur Wash 1^{a,d} 2 x 1500 mL/pack Phosphate-buffered saline; sodium azide (< 0.1%); surfactant | Unopened at 2–25° C | Until expiration date on product |
| | Onboard | 1 month |
| ADVIA Centaur Wash 1^{a,d} 2 x 2500 mL/pack Phosphate-buffered saline; sodium azide (< 0.1%); surfactant | Unopened at 2–25° C | Until expiration date on product |
| | Onboard | 1 month |

^a Store in an upright position.

^b Prevent exposure to sunlight and heat.

^c Refer to Optional Materials.

^d Refer to *Materials Required but Not Provided*.

Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

Safety data sheets (SDS) available on [siemens-healthineers.com](https://www.siemens-healthineers.com).



H317
H319
H412
P261, P280,
P273,
P302+P352,
P305+P351+
P338
P333+P313,
P362+P364,
P501

Warning!

May cause an allergic skin reaction.
Causes serious eye irritation.
Harmful to aquatic life with long lasting effects.
Avoid breathing vapor. Wear protective gloves/protective clothing/eye protection/face protection. Avoid release to the environment. IF ON SKIN: Wash with plenty of soap and water. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse. Dispose of contents and container in accordance with all local, regional, and national regulations.
Contains: reaction mass of 5-chloro-2-methyl-2H-isothiazol-3-one and 2-methyl-2H-isothiazol-3-one (3:1) (ADVIA Centaur SCC Solid Phase and Lite Reagent).

CAUTION

This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Contains sodium azide as a preservative. Sodium azide can react with copper or lead plumbing to form explosive metal azides. On disposal, flush reagents with a large volume of water to prevent buildup of azides. Disposal into drain systems must be in compliance with prevailing regulatory requirements.

Storage and Stability

Store all reagents in an upright position, away from light and heat. Do not use products beyond the expiration date printed on the product labeling.

For information about product storage and stability, refer to *Reagents*.

Onboard Stability

Discard products at the end of the onboard stability interval. Do not use products beyond the expiration date printed on the product labeling.

For information about product onboard stability, refer to *Reagents*.

Specimen Collection and Handling

Serum and plasma (EDTA) are the recommended specimen types for this assay.

The handling and storage information provided here is based on data or references maintained by the manufacturer. It is the responsibility of the individual laboratory to use all available references and/or its own studies when establishing alternate stability criteria to meet specific needs.

Collecting the Specimen

- Wear proper protective gear, including gloves and a face mask, to ensure SCC will not be shed into the specimen.
- Observe universal precautions when collecting specimens. Handle all specimens as if they are capable of transmitting disease.¹⁴
- Follow recommended procedures for collection of diagnostic blood specimens by venipuncture.¹⁵
- Follow the instructions provided with your specimen collection device for use and processing.¹⁶
- Allow blood specimens to clot completely before centrifugation.¹⁷
- Keep tubes capped at all times.¹⁷

Storing the Specimen

- Uncentrifuged specimens are stable under the following conditions:
 - Serum and tripotassium EDTA plasma in the primary collection device are stable for up to 3 days at 2–8° C.
 - Dipotassium EDTA plasma is unstable if left unprocessed. Do not delay centrifugation.
- After centrifugation, specimens left in the primary collection tube in contact with the cells or on the clot are stable under the following conditions:
 - Serum samples are stable for up to 4 days at room temperature, and for up to 4 days at 2–8° C.
 - Plasma samples are stable for up to 1 day at room temperature, and for up to 1 day at 2–8° C.
- Separated samples are stable under the following conditions:
 - Red Top serum samples are stable for up to 7 days at room temperature, and for up to 1 month at 2–8° C.
 - SST serum samples are stable for up to 7 days at room temperature, and for up to 1 day at 2–8° C.
 - Plasma samples are stable for up to 7 days at room temperature, and for up to 1 month at 2–8° C.
 - All samples are stable at $\leq -10^{\circ}$ C for up to 1 month. Avoid more than 1 freeze-thaw cycle. Do not store in a frost-free freezer. Thoroughly mix thawed samples and centrifuge them before using.

Transporting the Specimen

Package and label specimens for shipment in compliance with applicable federal and international regulations covering the transport of clinical specimens and etiological agents.

If during shipment, specimens may be subjected to temperatures $> 25^{\circ}$ C, then ship specimens frozen.

Preparing the Samples

This assay requires 50 µL of sample for a single determination. This volume does not include the unusable volume in the sample container or the additional volume required when performing duplicates or other tests on the same sample. For a complete list of appropriate sample containers and information about determining the minimum required volume, refer to the system online help.

Do not use samples with apparent contamination.

Before placing samples on the system, ensure that samples are free of:

- Bubbles or foam.
- Fibrin or other particulate matter.

Remove particulates by centrifugation according to CLSI guidance and the collection device manufacturer's recommendations.¹⁷

Procedure

Materials Provided

The following materials are provided:

| REF | Contents | Number of Tests |
|----------|--------------------------------------------------------------------------------------------|-----------------|
| 11354594 | 1 ReadyPack primary reagent pack containing ADVIA Centaur SCC Lite Reagent and Solid Phase | 100 |
| | ADVIA Centaur SCC master curve and test definition | MC TDEF |

Materials Required but Not Provided

The following materials are required to perform this assay, but are not provided:

| REF | Description |
|----------------------|------------------------------------------------|
| | ADVIA Centaur CP System ^a |
| 11553960 | ADVIA Centaur SCC CAL (calibrator) CAL LOT VAL |
| 01137199 (112351) | ADVIA Centaur Wash 1 (wash) WASH 1 |
| 03773025 | ADVIA Centaur Wash 1 (wash) WASH 1 |

^a Additional system fluids are required to operate the system: ADVIA Centaur Acid Reagent, ADVIA Centaur Base Reagent, and ADVIA Centaur Cleaning Solution.

Optional Materials

The following materials may be used to perform this assay, but are not provided:

| REF | Description | |
|----------|------------------------------------------|-----------------|
| 11538186 | ADVIA Centaur TM QC (quality controls) | CONTROL LOT VAL |
| 10492364 | ADVIA Centaur Multi-Diluent 13 (diluent) | DIL |

Assay Procedure

The system automatically performs the following steps:

1. Dispenses 50 µL of sample into a cuvette.
2. Dispenses 175 µL of Solid Phase, then incubates for 6 minutes at 37° C.
3. Dispenses 100 µL of Lite Reagent, then incubates for 6 minutes at 37° C.
4. Performs a wash sequence using ADVIA Centaur Wash 1.
5. Dispenses 300 µL each of ADVIA Centaur Acid Reagent and ADVIA Centaur Base Reagent to initiate the chemiluminescent reaction.
6. Reports results.

Preparing the Reagents

All reagents are liquid and ready to use. Before loading the packs onto the system, reagents require mixing. For information about mixing the reagents, refer to the system online help.

Preparing the System

Ensure that sufficient materials are loaded on the system. Refer to *Materials Provided*, *Materials Required but Not Provided*, and *Optional Materials* for guidance about required reagents.

For information about loading products, refer to the system online help.

Master Curve Definition

Before initiating calibration on each new lot of reagent, enter the assay master curve values by scanning the master curve card. For information about defining the master curve, refer to the system online help.

Performing Calibration

For calibration of the ADVIA Centaur SCC assay, use ADVIA Centaur SCC CAL. Use the calibrators in accordance with the calibrator instructions for use.

Calibration Frequency

Perform a calibration if one or more of the following conditions exist:

- At the end of the 30-day calibration interval.
- When changing lot numbers of primary reagent packs.

- When indicated by quality control results.
- After major maintenance or service, if indicated by quality control results.

Follow government regulations or accreditation requirements for calibration frequency. Individual laboratory quality control programs and procedures may require more frequent calibration.

Performing Quality Control

For quality control of the ADVIA Centaur SCC assay, use the ADVIA Centaur TM QC or an equivalent product at least once during each day that samples are analyzed. Use the quality control material in accordance with the quality control instructions for use. For the assigned values, refer to the quality control assigned value sheet provided.

Additional quality control material can be used at the discretion of the laboratory. Use the quality control material in accordance with the quality control instructions for use.

In addition, perform quality control:

- Following a valid calibration
- With use of a new lot of reagent
- When troubleshooting test results that do not match clinical conditions or symptoms

Follow government regulations or accreditation requirements for quality control frequency. Individual laboratory quality control programs and procedures may require more frequent quality control testing.

Acceptable performance is achieved when the analyte values obtained are within the expected control interval for the system, as indicated by the manufacturer of the control material or within the interval determined by an internal laboratory quality control procedure.

Follow your laboratory's quality control procedures if the results obtained do not fall within the acceptable limits. For information about entering quality control definitions, refer to the system online help.

Taking Corrective Action

If the quality control results do not fall within the expected control interval, do not report results. Perform corrective actions in accordance with established laboratory protocol. For suggested protocol, refer to the system online help.

Results

Calculation of Results

The system determines the result using the calculation procedure described in the system online help. The system reports results in ng/mL (common units) or µg/L (SI units), depending on the units defined when setting up the assay.

Conversion formula: 1 ng/mL = 1 µg/L

For information about results outside the specified measuring interval, refer to *Measuring Interval*.

Dilutions

The measuring interval is 0.60–70.00 ng/mL (0.60–70.00 µg/L). For information about dilution options, refer to the system online help.

Dilute and retest samples with SCC levels > 70.00 ng/mL (70.00 µg/L) to obtain accurate results.

For automated dilutions, perform the following activities:

- Load ADVIA Centaur Multi-Diluent 13.
- Ensure that sufficient sample volume is available. Refer to the table below.
- Select the appropriate dilution factor.

For automatic dilutions enter a Dilution Point ≤ 70 ng/mL (70 µg/L).

| Sample | Dilution | Sample Volume (µL) |
|------------------|----------|--------------------|
| Serum and plasma | 1:10 | 50 |

Interpretation of Results

Results of this assay should always be interpreted in conjunction with the patient's medical history, clinical presentation, and other findings.

Limitations

The following information pertains to limitations of the assay:

- Do not perform manual dilution.
- Dipotassium EDTA plasma is unstable if left unprocessed. Do not delay centrifugation.
- Do not use Lithium Heparin plasma tubes.
- Results obtained with the assay may not be used interchangeably with values obtained with different manufacturers' SCC methods.
- Patient samples may contain heterophilic antibodies that could react in immunoassays and cause falsely elevated or depressed results. This assay is designed to minimize interference from heterophilic antibodies.^{18,19} Additional information may be required for diagnosis.

Expected Values

Expected values were established using the ADVIA Centaur XPT system and confirmed by assay comparison. Refer to *Assay Comparison*.

A reference interval for healthy adults (smokers and nonsmokers) was established in accordance with CLSI Document EP28-A3c²⁰ on the ADVIA Centaur XPT system.

Samples tested were a combination of commercially available samples (N=14) and samples collected prospectively from apparently healthy subjects (N=226), males and females ≥ 22

years old without any chronic medical conditions or on chronic-use medications, women using hormones for menopause or on birth control were eligible. The reference interval was determined by calculating the 95th and 97.5th percentiles of the distribution of values.

| Group | N ^a | Median | Reference Interval |
|----------------------------|----------------|---------------|--------------------------------------------------------------------|
| | | ng/mL or µg/L | Percentile (CI) |
| Adults (22–87 years) | 240 | 1.24 | 2.86 (2.61 – 3.13) ^b 3.35 (3.01 – 3.71) ^c |
| Smoker (22–73 years) | 120 | 1.26 | 3.06 (2.66 – 3.47) ^b 3.61 (3.11 – 4.20) ^c |
| Nonsmoker (22–87 years) | 120 | 1.21 | 2.71 (2.39 – 3.06) ^b 3.17 (2.70 – 3.72) ^c |

^a Number of samples tested.

^b 95th percentile

^c 97.5th percentile

The distribution in percentage (%) of SCC assay values in benign and malignant cohorts was determined using 829 serum samples obtained from 9 U.S.-based clinical centers and commercially available sources using the ADVIA Centaur SCC assay.

| | | | Percentage | Mean | SD | 25 th | Median | 75 th |
|------------------------------|----------------|------|------------------|------------------|------------------|------------------|------------------|------------------|
| | | | ≥ 2.86 ng/mL | | | Percentile | | percentile |
| Group | N ^a | (%) | ng/mL or µg/L | ng/mL or µg/L | ng/mL or µg/L | ng/mL or µg/L | ng/mL or µg/L | ng/mL or µg/L |
| Benign Conditions | | | | | | | | |
| Lung Diseases | 75 | 18.7 | 2.27 | 2.64 | 1.05 | 1.43 | 2.43 | |
| Breast Diseases | 37 | 0.0 | 1.01 | 0.36 | 0.74 | 0.92 | 1.31 | |
| Liver Diseases | 38 | 2.6 | 1.39 | 0.61 | 0.83 | 1.34 | 1.76 | |
| Renal Diseases | 39 | 30.8 | 3.72 | 3.96 | 1.23 | 2.18 | 4.39 | |
| Benign Prostatic Hyperplasia | 40 | 17.6 | 3.06 | 4.70 | 1.24 | 1.81 | 2.67 | |
| Cancer | | | | | | | | |
| Treatment-naïve NSCLC | 120 | 2.5 | 1.34 | 0.95 | 0.76 | 1.17 | 1.69 | |
| Bladder | 40 | 20.0 | 2.87 | 3.59 | 1.12 | 1.77 | 2.53 | |
| Breast | 48 | 2.1 | 1.23 | 0.77 | 0.77 | 0.96 | 1.60 | |
| Cervical | 40 | 35.0 | 8.21 | 15.34 | 1.06 | 1.58 | 5.28 | |
| Colorectal | 40 | 15.0 | 2.75 | 6.74 | 0.98 | 1.42 | 2.17 | |
| Esophageal SCC | 35 | 28.6 | 2.78 | 1.88 | 1.45 | 2.07 | 3.43 | |
| Head & Neck | 40 | 15.0 | 4.22 | 12.56 | 1.12 | 1.53 | 2.22 | |

| Group | N ^a | Percentage ≥ 2.86 ng/mL | Mean | SD | 25 th Percentile | Median | 75 th percentile |
|----------------|----------------|----------------------------|------------------|------------------|--------------------------------|------------------|--------------------------------|
| | | (%) | ng/mL or µg/L | ng/mL or µg/L | ng/mL or µg/L | ng/mL or µg/L | ng/mL or µg/L |
| Neuroendocrine | 37 | 10.8 | 1.43 | 0.78 | 0.97 | 1.19 | 1.59 |
| Ovarian | 40 | 5.0 | 1.39 | 1.46 | 0.71 | 0.91 | 1.55 |
| Prostate | 40 | 12.5 | 1.72 | 1.19 | 0.92 | 1.39 | 1.97 |
| Renal | 40 | 20.0 | 2.59 | 3.80 | 1.23 | 1.90 | 2.56 |
| Stomach | 40 | 5.0 | 1.50 | 1.04 | 0.91 | 1.22 | 1.68 |
| Testicular | 40 | 7.5 | 1.54 | 0.87 | 0.86 | 1.37 | 1.94 |

^a Number of samples tested.

As with all *in vitro* diagnostic assays, each laboratory should determine its own reference interval for the diagnostic evaluation of patient results.²⁰ Consider these values as guidance only.

Performance Characteristics

Measuring Interval

0.60–70.00 ng/mL (0.60–70.00 µg/L)

The lower limit of the measuring interval is defined by the limit of quantitation (LoQ). Report results below the measuring interval as < 0.60 ng/mL (0.60 µg/L).

When sample results exceed the measuring interval, refer to *Dilutions*.

Detection Capability

Detection capability was determined in accordance with CLSI Document EP17-A2.²¹

| Method | Result |
|-----------------------------|---------------|
| | ng/mL or µg/L |
| Limit of Blank (LoB) | 0.18 |
| Limit of Detection (LoD) | 0.34 |
| Limit of Quantitation (LoQ) | 0.34 |

These are representative data and the assay results obtained at individual laboratories may vary from the data presented.

The LoB corresponds to the highest measurement result likely to be observed for a blank sample with a probability of 95%. The assay was designed to have an LoB of < 0.60 ng/mL (0.60 µg/L).

The LoD corresponds to the lowest concentration of SCC Ag that can be detected with a probability of 95%. The assay was designed to have an LoD of ≤ 0.60 ng/mL (0.60 µg/L).

The LoQ corresponds to the lowest amount of SCC Ag in a sample at which the within-laboratory precision (%CV) is $\leq 20\%$. The assay was designed to have an LoQ of ≤ 0.60 ng/mL (0.60 µg/L).

Precision

Precision was determined in accordance with CLSI Document EP05-A3.²² Samples were assayed in duplicate in 2 runs per day for 20 days.

The following results are representative of the performance of the assay:

| Specimen Type | N ^a | Repeatability | | | Within-Laboratory Precision | |
|---------------|----------------|---------------|-----------------|------------------|-----------------------------|------------------|
| | | Mean | SD ^b | CV ^c | SD | CV |
| | | ng/mL or µg/L | ng/mL or µg/L | (%) | ng/mL or µg/L | (%) |
| Serum A | 80 | 0.90 | 0.042 | N/A ^d | 0.046 | N/A ^d |
| Serum B | 80 | 3.20 | 0.087 | 2.7 | 0.198 | 3.1 |
| Serum C | 80 | 8.36 | 0.169 | 2.0 | 0.254 | 3.0 |
| Serum D | 80 | 38.30 | 0.547 | 1.4 | 1.005 | 2.6 |
| Serum E | 80 | 66.00 | 1.137 | 1.7 | 1.332 | 2.0 |

^a Number of measurements.

^b Standard deviation.

^c Coefficient of variation.

^d Not applicable.

The assay was designed to have the following precision.

| Concentration Interval | Precision | |
|------------------------|------------------------------------------------|------------------------------------------------|
| ng/mL or µg/L | Repeatability (Within-Run) | Within-Laboratory (Total Precision) |
| ≥ 3.00 | $\leq 7.5\%$ CV | $\leq 7.5\%$ CV |
| ≥ 0.41 - < 3.00 | ≤ 0.23 ng/mL SD (≤ 0.23 µg/L SD) | ≤ 0.23 ng/mL SD (≤ 0.23 µg/L SD) |

Assay Comparison

Assay comparison was determined with the weighted Deming regression model in accordance with CLSI Document EP09c-ed3.²³

Agreement of the assays may vary depending on the study design, comparative assay, and population tested.

| Specimen Type | Comparative Assay (x) | Regression Equation | Sample Interval | | |
|---------------|------------------------------------------------------------|-------------------------------------------------------------------------|-----------------|----------------|----------------|
| | | | ng/mL or µg/L | N ^a | r ^b |
| Serum | ADVIA Centaur SCC assay using the ADVIA Centaur XPT system | $y = 1.05x - 0.02 \text{ ng/mL}$ ($y = 1.05x - 0.02 \text{ µg/L}$) | 0.55–54.61 | 119 | 0.997 |

^a Number of samples tested.

^b Correlation coefficient.

The assay is designed to have a correlation coefficient of ≥ 0.97 , a slope of 1.00 ± 0.10 , and an intercept within $\pm 0.60 \text{ ng/mL}$ (0.60 µg/L).

Specimen Equivalency

Specimen equivalency was determined with the weighted Deming regression model using the ADVIA Centaur XPT system in accordance with CLSI Document EP09c-ed3.²³

Agreement of the specimen types may vary depending on the study design and population tested.

| Tube (y) vs. Serum (x) | Regression Equation | Sample Interval | | |
|----------------------------|-------------------------------------------------------------------------|-----------------|----------------|----------------|
| | | ng/mL or µg/L | N ^a | r ^b |
| Dipotassium EDTA plasma | $y = 1.02x + 0.08 \text{ ng/mL}$ ($y = 1.02x + 0.08 \text{ µg/L}$) | 0.78–62.62 | 50 | 0.997 |
| Tripotassium EDTA plasma | $y = 0.98x + 0.08 \text{ ng/mL}$ ($y = 0.98x + 0.08 \text{ µg/L}$) | 0.64–69.16 | 50 | 0.998 |
| Serum separator tube (SST) | $y = 0.99x - 0.07 \text{ ng/mL}$ ($y = 0.99x - 0.07 \text{ µg/L}$) | 0.68–61.65 | 50 | 0.996 |

^a Number of samples tested.

^b Correlation coefficient.

The assay is designed to have a correlation coefficient of ≥ 0.97 , a slope of 1.0 ± 0.10 , and an intercept within $\pm 0.60 \text{ ng/mL}$ (0.60 µg/L).

Interferences

Hemolysis, Icterus, Lipemia (HIL)

Interference testing was performed using the ADVIA Centaur XPT system in accordance with CLSI Document EP07-ed3.²⁴

The following substances do not interfere with the assay when present in serum at the concentrations indicated. Bias due to these substances does not exceed 10% at analyte concentrations of approximately 2.50 ng/mL (2.50 µg/L) and 50.0 ng/mL (50.0 µg/L).

| Substance | Substance Test Concentration |
|-------------------------|------------------------------|
| Hemoglobin | 1012 mg/dL |
| Bilirubin, conjugated | 60 mg/dL |
| Bilirubin, unconjugated | 67 mg/dL |
| Lipemia (Triglycerides) | 1516 mg/dL |
| Lipemia (Intralipid) | 1500 mg/dL |

Other Substances

Interference testing was performed using the ADVIA Centaur XPT system in accordance with CLSI Document EP07-ed3²⁴ and EP37-ed1.²⁵

The following substances do not interfere with the assay when present in serum at the concentrations indicated. Bias due to these substances does not exceed 10% at analyte concentrations of approximately 2.50 ng/mL (2.50 µg/L) and 50.0 ng/mL (50.0 µg/L).

| Substance | Substance Test Concentration |
|----------------------------------------|------------------------------|
| Total Protein | 15 g/dL |
| Immunoglobulin G | 20 g/L |
| Immunoglobulin A | 5 g/L |
| Immunoglobulin M | 10 g/L |
| Rheumatoid Factor | 1003 IU/mL |
| Human Anti-Mouse Antibodies (HAMA IgG) | 1000 µg/L |
| 5-Fluorouracil | 900 µg/mL |
| Acetaminophen | 200 µg/mL |
| Acetylcysteine | 553 µg/mL |
| Acetylsalicylic acid | 1000 µg/mL |
| Afatinib | 0.024 mg/mL |
| Ampicillin-Na | 1000 µg/mL |
| Aprepitant | 0.075 mg/mL |
| Ascorbic acid | 300 µg/mL |
| Atezolizumab | 1.008 mg/mL |
| Bevacizumab | 700 µg/mL |
| Biotin | 3500 ng/mL |
| Carboplatin | 600 µg/mL |
| Cefoxitin | 2500 µg/mL |
| Cetuximab | 600 µg/mL |
| Cisplatin | 180 µg/mL |

| Substance | Substance Test Concentration |
|-------------------------|------------------------------|
| Cyclophosphamide | 500 µg/mL |
| Cyclosporine | 5 µg/mL |
| Dexamethasone | 20 µg/mL |
| Docetaxel | 112.5 µg/mL |
| Doxorubicin | 120 µg/mL |
| Doxycycline | 50 µg/mL |
| Epoetin Alfa | 0.378 µg/mL |
| Erlotinib | 150 µg/mL |
| Etoposide | 300 µg/mL |
| Fosaprepitant | 0.09 mg/mL |
| Gefitinib | 250 µg/mL |
| Gemcitabine | 1500 µg/mL |
| Heparin | 5000 U/L |
| Ibuprofen | 500 µg/mL |
| Ifosfamide | 7200 µg/mL |
| Levodopa | 20 µg/mL |
| Methotrexate | 150 µg/mL |
| Methyldopa +1.5 | 20 µg/mL |
| Metoclopramide | 7.5 µg/mL |
| Metronidazole | 200 µg/mL |
| Nab-Paclitaxel | 0.06 mg/mL |
| Neupogen | 0.9 µg/mL |
| Nivolumab | 0.288 mg/mL |
| Ondansetron | 0.0342 mg/dL |
| Paclitaxel | 265 µg/mL |
| Palonosetron | 0.00015 mg/mL |
| Pegfilgrastim | 0.0036 mg/mL |
| Pembrolizumab | 0.24 mg/mL |
| Phenylbutazone | 400 µg/mL |
| Prochlorperazine | 0.345 mg/dL |
| Rifampicin | 60 µg/mL |
| Theophylline | 100 µg/mL |
| Topotecan hydrochloride | 2.25 µg/mL |
| Vincristine sulfate | 3 µg/mL |
| Vinorelbine tartrate | 53.1 µg/mL |

Cross-Reactivity

Cross-reactivity was determined using the ADVIA Centaur XPT system in accordance with CLSI Document EP07-ed3.²⁴

| Cross-reactant | Cross-reactant Concentration | SCC concentration | Cross-reactivity |
|---------------------|------------------------------|-------------------|------------------|
| | | (ng/mL or µg/L) | (%) |
| ProGRP ^a | 1900 ng/L | 2.50 50.00 | 8 21 |
| CA 125 | 950 U/mL | 2.50 50.00 | < 1 < 1 |
| NSE | 400 ng/mL | 2.50 50.00 | < 1 < 1 |
| CYFRA 21-1 | 100 ng/mL | 2.50 50.00 | < 1 1 |
| CEA | 70 µg/L | 2.50 50.00 | < 1 2 |

^a Due to the extremely low concentration, cross reactivity for ProGRP can appear unusually high. The Percent Difference for ProGRP ranged from 0% to 7%.

Linearity

Linearity testing was performed in accordance with CLSI Document EP06-ed2.²⁶

Linearity was evaluated using a sample that contained a high level of SCC, mixed in various proportions with a sample that contained a low level of SCC. The resulting sample mixtures (15 combinations) were assayed for SCC.

The ADVIA Centaur SCC assay is linear for the measuring interval of 0.60–70.00 ng/mL (0.60–70.00 µg/L).

Onboard Dilution Recovery

Serum samples were diluted onboard the ADVIA Centaur CP system with ADVIA Centaur Multi-Diluent 13.

The following results are representative of the performance of the assay:

| Sample | Dilution | Observed | Expected | Recovery |
|--------|----------|---------------|---------------|----------|
| | | ng/mL or µg/L | ng/mL or µg/L | (%) |
| 1 | 1:10 | 82.54 | 85 | 97 |
| 2 | 1:10 | 190.69 | 200 | 95 |
| 3 | 1:10 | 348.93 | 350 | 100 |
| 4 | 1:10 | 507.28 | 500 | 101 |
| 5 | 1:10 | 658.90 | 650 | 101 |
| Mean | | | | 99 |

High-Dose Hook Effect

High SCC concentrations can cause a paradoxical decrease in the RLUs (high-dose hook effect). In this assay, patient samples with SCC concentrations above the measuring interval and as high as 5,000 ng/mL (5,000 µg/L) will report >70.00 ng/mL (70.00 µg/L).

Standardization

The ADVIA Centaur SCC assay is traceable to an internal standard manufactured using highly purified material.

Currently no reference standard is available for this assay.

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

For customer support, contact your local technical support provider or distributor.

[siemens-healthineers.com](https://www.siemens-healthineers.com)

References









1. Kato, H, Miyauchi, F, Morioka, H, et al. Tumor antigen of human cervical squamous cell carcinoma: correlation of circulating levels with disease progress. *Cancer*. 1979;43(2):585-590.
2. Schneider, SS, C, Fish KE, et al. A serine proteinase inhibitor locus at 18q21. 3 contains a tandem duplication of the human squamous cell carcinoma antigen gene. *Proc Natl Acad Sci USA* 1995;92(8):3147-3151.



















3. Röijer, E, Kosinska U, Andersson I, et al. Rearrangement of squamous cell carcinoma antigen genes—detection of SCCA fusion transcripts. *Tumour Biol.* 2003;24(1):46-52.
4. Cataltepe, S, Gornstein ER, Schick C, et al. Co-expression of the squamous cell carcinoma antigens 1 and 2 in normal adult human tissues and squamous cell carcinomas. *J Histochem Cytochem.* 2000;48(1):113-122.
5. Yagi, H, Danno K, Maruguchi Y, et al. Significance of squamous cell carcinoma (SCC)-related antigens in cutaneous SCC: A preliminary report. *Arch Dermatol.* 1987;123(7):902-906.
6. Yoshimura, Y, Harada T, Oka M, et al. Squamous cell carcinoma antigen in the serum of oromaxillary cancer. *Int J Oral Maxillofac Surg.* 1988;17(1):49-53.
7. Petrelli, NJ, Shaw N, Bhargava A, et al. Squamous cell carcinoma antigen as a marker for squamous cell carcinoma of the anal canal. *J Clin Oncol.* 1988;6(5):782-785.
8. Yamanaka N, Himi T, Harabuchi Y, et al. Soluble immune complexes and squamous cell carcinoma-related antigens in patients with head and neck cancer. *Cancer.* 1988;62(9):1932-1938.
9. Mino N, Iio A, Hamamoto K. Availability of tumor-antigen 4 as a marker of squamous cell carcinoma of the lung and other organs. *Cancer.* 1988;62(4):730-734.
10. Eibling DE, Johnson JT, Wagner RL, et al. SCC-RIA in the diagnosis of squamous cell carcinoma of the head and neck. *Laryngoscope.* 1989;99(2):117-124.
11. Johnson J, et al. Radioimmunoassay for SCC antigen in the diagnosis of squamous cell carcinoma of the head and neck: A preliminary report. *SCC antigen in the management of squamous cell carcinoma. Princeton: Excerpta Medica.* 1987:112-123.
12. Bonfrer JMG, Duffy MJ, Radtke M, et al. Tumour Markers in Gynaecological Cancers - EGTm Recommendations. *Anticancer Research.* 1999;19(4A):2807-2810.
13. Fu J, Wang W, Wang Y, et al. The role of squamous cell carcinoma antigen (SCC Ag) in outcome prediction after concurrent chemoradiotherapy and treatment decisions for patients with cervical cancer. *Radiat Oncol.* 2019;14(1):146.
14. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
15. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
16. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
17. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
18. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem.* 1999;45(7):942-956.
19. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem.* 1992;38(9):1737-1742.


















20. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
21. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
22. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
23. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
24. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP07-ed3.
25. Clinical and Laboratory Standards Institute. *Supplemental Tables for Interference Testing in Clinical Chemistry; Approved Supplement—First Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP37.
26. Clinical and Laboratory Standards Institute. *Evaluation of Linearity of Quantitative Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2020. CLSI Document EP06-ed2.

Definition of Symbols

The following symbols may appear on the product labeling:

| Symbol | Symbol Title and Description |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Consult instructions for use |
|  Rev. 01 | Version of instructions for use |
|  siemens.com/healthcare | Internet URL address to access the electronic instructions for use |
|  siemens.com/document-library | |
| Rev. REVISION | Revision |
|  | Caution Consult instructions for use or accompanying documents for cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device. |
|  | Biological risks Potential biological risks are associated with the medical device. |
|  | Corrosive |
|  | Dangerous to environment |

| Symbol | Symbol Title and Description |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Irritant Oral, dermal, or inhalation hazard |
|  | Inhalation hazard Respiratory or internal health |
|  | Flammable Flammable to extremely flammable |
|  | Oxidizing |
|  | Explosive |
|  | Toxic |
|  | Compressed gas |
|  | Keep away from sunlight Prevent exposure to sunlight and heat. |
|  | Up Store in an upright position. |
|  | Temperature limit Upper and lower limits of temperature indicators are adjacent to the upper and lower horizontal lines. |
|  | Do not freeze. |
|  | Handheld barcode scanner |
|  | <i>In vitro</i> diagnostic medical device |
|  | Contains sufficient for <n> tests Total number of IVD tests the system can perform with the IVD kit reagents appears adjacent to the symbol. |
| RxOnly | Prescription device (US only) Applies only to United States-registered IVD assays. CAUTION: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional. |
|  | Mixing of substances Mix product before use. |
|  | Target |
|  | Interval |
|  | Legal Manufacturer |

| Symbol | Symbol Title and Description |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
|  | Authorized Representative in the European Community |
|  | Use-by date Use by the designated date. |
|  | Batch code |
|  | Catalog number |
|  | Recycle |
|  | Printed with soy ink |
|  | CE Mark |
|  | CE Mark with notified body ID number Notified body ID number can vary. |
| XXXX | |
| YYYY-MM-DD | Date format (year month day) |
|  | Variable hexadecimal number that ensures the Master Curve and Calibrator definition values entered are valid. |
|  | Master Curve Definition |
|  | Lot Details |
|  | Common Units |
|  | International System of Units |
|  | Material |
|  | Unique material identification number |
|  | Name of control |
|  | Type of control |

Legal Information

ADVIA Centaur and ReadyPack are trademarks of Siemens Healthineers.

All other trademarks and brands are the property of their respective owners.

© 2022 Siemens Healthineers. All rights reserved.

US Pats 9,575,062; 6,664,043; 8,778,624; 7,309,615; 7,785,904.



Siemens Healthcare Diagnostics Inc.

511 Benedict Avenue

Tarrytown, NY 10591 USA

Siemens Healthineers Headquarters:

Siemens Healthcare GmbH

Heenkestraße 127

91052 Erlangen

Germany

Phone: +49 9131 84-0

[siemens-healthineers.com](https://www.siemens-healthineers.com)