

Calibrator Q (CAL Q)

Current Revision and Date^a	Rev. 03, 2022-10
Product Name	Atellica IM Calibrator Q (CAL Q)
Abbreviated Product Name	Atellica IM CAL Q
	<div>2 x 2.0 mL low calibrator CAL L REF 10995517</div> <div>2 x 2.0 mL high calibrator CAL H (2-pack)</div> <div>Calibrator lot-specific value sheet CAL LOT VAL</div> <div>6 x 2.0 mL low calibrator CAL L REF 10995518</div> <div>6 x 2.0 mL high calibrator CAL H (6-pack)</div> <div>Calibrator lot-specific value sheet CAL LOT VAL</div>
Systems	Atellica IM Analyzer

^a A vertical bar in the page margin indicates technical content that differs from the previous version.



Intended Use

The Atellica® IM Calibrator Q (CAL Q) is for *in vitro* diagnostic use in calibrating the Atellica IM PSA assay using the Atellica® IM Analyzer.

Material Description

For the assigned values, refer to the calibrator lot-specific value sheet CAL LOT VAL provided. The assigned values are traceable to the standardization of the assay. For additional information, refer to the assay instructions for use.

Material Description	Storage	Stability ^a
Atellica IM CAL Q 2.0 mL/vial; lyophilized After reconstitution, low or high levels of PSA (human); goat serum; sodium azide (< 0.1%); preservatives	Lyophilized at 2–8°C	Until expiration date on product
	Reconstituted at 2–8°C	21 days
	Reconstituted at room temperature	8 hours

^a Refer to *Storage and Stability*.

Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

CAUTION

Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

Safety data sheets (SDS) available on [siemens-healthineers.com](https://www.siemens-healthineers.com).

The summary of safety and performance for this *in vitro* diagnostic medical device is available to the public in the European database on medical devices (EUDAMED) when this database is available and the information has been uploaded by the Notified Body. The web address of the EUDAMED public website is: <https://ec.europa.eu/tools/eudamed>.

H412	Harmful to aquatic life with long lasting effects.
P273, P501	Avoid release to the environment. Dispose of contents and container in accordance with all local, regional, and national regulations.
	Contains: sodium azide (Atellica IM CAL Q)

**Warning! Potential Biohazard**

Contains human source material.

No known test method can ensure that products derived from human source materials will not transmit infection. These materials should be handled using good laboratory practices and universal precautions.¹⁻³

CAUTION

This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

Contains sodium azide as a preservative. Sodium azide can react with copper or lead plumbing to form explosive metal azides. On disposal, flush reagents with a large volume of water to prevent buildup of azides. Disposal into drain systems must be in compliance with prevailing regulatory requirements.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.

Note For information about calibrator preparation, refer to *Preparing the Calibrators*.

Storage and Stability

Store calibrators in an upright position. Lyophilized calibrators are stable until the expiration date on the product when stored at 2–8°C. Reconstituted calibrators are stable for 21 days at 2–8°C. Reconstituted calibrators are stable for 8 hours at room temperature.

Do not use products beyond the expiration date printed on the product labeling.

Performing Calibration

Calibration Frequency

For information about calibration frequency, refer to the assay instructions for use.

Preparing the Calibrators

Prepare calibrators using the following steps:

1. Add 2.0 mL of special reagent water into each vial using a precision pipet. Replace cap.
Note For information about special reagent water requirements, refer to the online help.
2. Let the vials stand for 15–20 minutes at room temperature to allow the lyophilized material to dissolve.
3. Gently mix and invert the vials to ensure homogeneity of the material.

Note Use calibrators within the stability limits specified in *Storage and Stability* and discard any remaining material.

Calibration Procedure

The required sample volume for testing depends on several factors. For information about sample volume requirements, refer to the online help.

Use the following lot-specific materials to perform calibration:

- For the master curve and assay test definitions, refer to the lot-specific master curve and test definition sheet **MC TDEF** provided with the assay reagents.
- For the calibrator definitions, refer to the lot-specific value sheet **CAL LOT VAL** provided with the calibrator materials.
- Generate lot-specific barcode labels to use with the calibrator samples.

For instructions about how to perform the calibration procedure, refer to the online help.

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

For customer support, contact your local technical support provider or distributor.



































siemens-healthineers.com

References

1. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
2. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
3. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.

Definition of Symbols

The following symbols may appear on the product labeling:

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Manufacturer	5.1.1 ^a		Authorized representative in the European Community	5.1.2 ^a
	Use-by date	5.1.4 ^a		Authorized representative in Switzerland	Proprietary
	Catalog number	5.1.6 ^a		Batch code	5.1.5 ^a
	Consult Instructions for Use	5.4.3 ^a		Contains sufficient for <n> tests	5.5.5 ^a
	Internet URL address to access the electronic instructions for use	Proprietary		Version of Instructions for Use	Proprietary
	<i>In vitro</i> diagnostic medical device	5.5.1 ^a		Revision	Proprietary
	Prescription device (US only)	FDA ^b		Unique Device Identifier	5.7.10 ^c
	CE Marking with Notified Body	EU IVDR ^d		CE Marking	EU IVDR ^d
	Temperature limit	5.3.7 ^a		Keep away from sunlight	5.3.2 ^a
	Upper limit of temperature	5.3.6 ^a		Lower limit of temperature	5.3.5 ^a
	Do not re-use	5.4.2 ^a		Do not freeze	Proprietary
	Recycle	1135 ^e		This way up	0623 ^e
	Biological risks	5.4.1 ^a		Caution	5.4.4 ^a
	Common Units	Proprietary		International System of Units	Proprietary
	Date format (year-month-day)	N/A		Date format (year-month)	N/A
	Document face up ^f	1952 ^e		Handheld barcode scanner	Proprietary
	Target	Proprietary		Mixing of substances	5657 ^g

Symbol	Symbol Title	Source	Symbol	Symbol Title	Source
	Variable hexadecimal number that ensures the Master Curve and Calibrator definition values entered are valid.	Proprietary		Interval	Proprietary
	Unique material identification number	Proprietary		Material	Proprietary
	Type of control	Proprietary		Name of control	Proprietary
	Quality control lot value	Proprietary		Calibrator lot value	Proprietary

- ^a International Standard Organization (ISO). ISO 15223-1 Medical Devices- Symbols to be used with medical device labels, labelling and information to be supplied.
- ^b Federal Register. Vol. 81, No 115. Wednesday, June 15, 2016. Rules and Regulations: 38911.
- ^c ISO 15223-1:2020-04
- ^d IVDR REGULATION (EU) 2017/746
- ^e International Standard Organization (ISO). ISO 7000 Graphical symbols for use on equipment.
- ^f Indicates Assay-eNote
- ^g International Electrotechnical Commission (IEC). IEC 60417-1 Graphical symbols for use on equipment – Part 1: Overview and Application

Legal Information

Atellica is a trademark of Siemens Healthineers.

All other trademarks and brands are the property of their respective owners.

© 2017–2022 Siemens Healthineers. All rights reserved.



Siemens Healthcare Diagnostics Inc.
511 Benedict Avenue
Tarrytown, NY 10591 USA

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestraße 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com

