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- Kaplan MM. Thyroid function testing in patients with thyroid and non-thyroid diseases. Mono: Thyroid Testing. Emeryville, CA: Chiron Diagnostics Corporation; 1996.
- Fisher DA. Physiological variations in thyroid hormones: physiological and pathophysiological considerations. Clin Chem. 1996;42(1):135-139.
- Chen IW, Sperling MI. Thyroxine. In: Kaplan LA, Pesce AJ, eds. Clinical Chemistry: Theory, Analysis, and Correlation. 2nd ed. St. Louis, MO: CV Mosby; 1989:952-956.
- Watts NB, Keffer JH. Practical Endocrine Diagnosis. 3rd ed. Philadelphia, PA: Lea and Febiger; 1982:1-27, 77-96.
- Chattoraj SC, Watts NB. Endocrinology. In: Tietz NW, ed. Fundamentals of Clinical Chemistry. 3rd ed. Philadelphia, PA: WB Saunders; 1986:550-551.
- Clinical and Laboratory Standards Institute. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Kricka LJ. Human anti-animal antibody interferences in

immunological assays. Clin Chem. 1999;45(7):942-956.

- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. Clin Chem. 1992;38(9):1737-1742.
- Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c (formerly C28-A3c).
- Horn PS, Pesce AJ. Reference Intervals: A User's Guide, Washington, DC: AACC Press; 2005.
- Clinical and Laboratory Standards Institute. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2000. CLSI Document C28-A2.
- Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
- Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.

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Siemens Healthcare Diagnostics Inc.