

Washing Solution for Coagulation Analyzers

WASH COAG

Revision bar indicates update to previous version.

Intended Use

Washing Solution is used for cleaning the Siemens Healthineers coagulation analyzers.

Reagents

Note: **WASH COAG** can be used on automated coagulation analyzers. Siemens Healthineers provides Reference Guides (Application Sheets) for several coagulation analyzers. The Reference Guides (Application Sheets) contain analyzer/assay specific handling and performance information which may differ from that provided in these Instructions for Use. In this case, the information contained in the Reference Guides (Application Sheets) supersedes the information in these Instructions for Use. Please also consult the instruction manual of the instrument manufacturer!

Reagent	Description	Storage	Stability
Washing Solution for Coagulation Analyzers WASH COAG	Ready to use liquid containing: <ul style="list-style-type: none"> Alcohols, C11-15-secondary, ethoxylated (Tergitol) (0.954 % [w/w]) <p>pH ≤ 2.0</p>	2–25 °C May be used up to the expiry date indicated on the label if stored unopened.	15–25 °C: once opened, 5 days

Warnings and Precautions

For *in-vitro* diagnostic use only.

For laboratory professional use.

Safety data sheets (MSDS/SDS) available on [siemens-healthineers.com/sds](https://www.siemens-healthineers.com/sds).

CAUTION!

Federal (USA) law restricts this device to sale by or on the order of licensed healthcare professionals.



Warning! **WASH COAG**

Hazardous ingredient: Hydrochloric acid (0.414 % [w/w]), Alcohols, C11-15-secondary, ethoxylated (Tergitol) (0.954 % [w/w]).

H290: May be corrosive to metals. **H319:** Causes serious eye irritation.

P234: Keep only in original container. **P280:** Wear protective gloves/protective clothing/eye protection/face protection. **P390:** Absorb spillage to prevent material damage. **P337 + P313:** If eye irritation persists: Get medical advice/attention.

Follow appropriate laboratory practices when handling solutions containing hydrochloric acid and Tergitol.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with all government requirements.

Preparing Reagents

Ready-for-use.

Procedure

Materials Provided

REF	Contents		
OWZC39	Washing Solution for Coagulation Analyzers WASH COAG	10 ×	15 mL

Materials Required but not Provided

























Item	Description
Coagulation analyzers, such as:	<ul style="list-style-type: none">BCS® XP System

Technical Assistance

For customer support, contact your local technical support provider or distributor.
siemens-healthineers.com

Definition of Symbols

The following symbols may appear on the product labeling:

	Do not reuse		Use By
	Batch Code		Catalogue Number
	Caution		Manufacturer
	Authorized representative in the European Community		Contains sufficient for <n> tests
	Biological Risks		<i>In Vitro</i> Diagnostic Medical Device
	Temperature Limitation		Consult instruction for Use
	Non-sterile		CE marking of conformity
	CE marking of conformity with notified body ID number. Notified body ID number can vary.		Contents
	Reconstitution volume		Level
	Keep away from sunlight and heat		Warning
	Danger		Prescription device (US only)
	Device Identification (UDI) barcode		REACH Authorization Number

Legal Information

BCS is a trademark of Siemens Healthineers.

All other trademarks are the property of their respective owners.

© Siemens Healthineers, 2008–2022. All rights reserved.

Siemens Healthineers Headquarters

Siemens Healthcare GmbH
Henkestraße 127
91052 Erlangen
Germany
Phone: +49 9131 84-0
siemens-healthineers.com



Siemens Healthcare Diagnostics Products GmbH

Emil-von-Behring-Str. 76
35041 Marburg
Germany
siemens-healthineers.com