IMMULITE®/IMMULITE® 1000

Chemiluminescent Substrate Module

English

Chemiluminescent Substrate Module

Intended Use: For use with the IMMULITE and IMMULITE 1000 Systems. Follow instructions in the IMMULITE or IMMULITE 1000 Operator's Manual.

Catalog Number: **LSUBX**

Summary and Explanation

Chemiluminescent Substrate is a dioxetane-based compound that is intended to be used in all IMMULITE assay reactions to generate light emission for qualitative or quantitative specific-analyte determinations.

Materials Supplied

Two bottles (LSUB5), 105 mL each (approximately 500 tests), of ready-to-use substrate, consisting of a phosphate ester of adamantyl dioxetane¹ in an AMP² buffer, with enhancer.³ Store at 2 – 8°C before opening. Stable on board after opening at 15 – 28°C for 30 days, or until expiration date.

Allow substrate to equilibrate at room temperature (about 20 minutes) prior to opening bottle and dispensing.

1. LUMIGEN® PPD: 4-methoxy-4-(3-phosphatephenyl)-spiro- (1,2-dioxetane-3,2'-adamantane). LUMIGEN® is a registered trademark of LUMIGEN, Inc., Southfield, MI. Patent numbers: 4,857,652, 4,983,779, 4,962,192, 5,386,017, and 5,707,559 (US); 254051B1 (Europe); 5-45590 (Japan).
2. AMP: 2-amino-2-methyl-1-propanol.
3. Polyvinylbenzyltributylphosphonium chloride. Patent numbers: 5,393,469 and 5,439,617 (US); 2648423 (Japan).

Risk and Safety

Safety data sheets (MSDS/SDS) available on siemens-healthineers.com

Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

Appearance

Precipitates or turbidity may indicate deterioration. Do not use if the normal appearance (clear and colorless) has been altered.

Light Sensitivity

Protect the Chemiluminescent Substrate from direct sunlight.

Avoid Contamination

Do not touch the spike with bare or gloved hands when changing bottles. Do not remove rubber stopper. Do not introduce syringe needles or other foreign objects into the bottle. Do not pool partial bottles of Substrate. If the bottle is removed from the IMMULITE Analyzer for any length of time, replace the protective plastic cover on the spike.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

In the United States, contact Siemens Healthcare Diagnostics Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

siemens-healthineers.com

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Shading indicates technical content that differs from the previous version.

IMMULITE® is a trademark of Siemens Healthineers.

All other trademarks and brands are the property of their respective owners.

© 2009–2022 Siemens Healthineers.
All rights reserved.

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