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siemens-healthineers.com

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1. Ruoslahti E, Seppala M. Studies of carcino-fetal proteins: physical and chemical properties of human alpha-fetoprotein. *Int J Cancer*. 1971;7(2):218–225.
2. Bergstrand CG, Czar B. Demonstration of a new protein fraction in serum from the human fetus. *Scand J Clin Lab Invest*. 1956;8(2):174.
3. Ruoslahti E, Engvall E, Kessler MJ. Chemical properties of alpha-fetoprotein. In: Herberman RB, McIntire KR, eds. *Immunodiagnosis of Cancer*. New York, NY: Marcel Dekker; 1979:101–117.
4. Morinaga T, Sakai M, Wegmann GT, Tamaoki T. Primary structures of human alpha-fetoprotein and its mRNA. *Proc Natl Acad Sci USA*. 1983;80(15):4604–4608.
5. Gitlin D, Perricelli A, Gitlin GM. Synthesis of α -fetoprotein by liver, yolk sac, and gastrointestinal tract of the human conceptus. *Cancer Res*. 1972;32(5):979–982.
6. Harris R, Jennison RF, Barson AJ, et al. Comparison of amniotic-fluid and maternal serum alpha-fetoprotein levels in the early antenatal diagnosis of spina bifida and anencephaly. *Lancet*. 1974;1(7855):428–429.
7. Brock DJ, Bolton AE, Monaghan JM. Prenatal diagnosis of anencephaly through maternal serum-alpha-fetoprotein measurement. *Lancet*. 1973;2(7835):923–924.
8. Wald NJ, Brock DJ, Bonnar J. Prenatal diagnosis of spina bifida and anencephaly by maternal serum-alpha-fetoprotein measurement. A controlled study. *Lancet*. 1974;1(7861):765–767.
9. Wald NJ, Cuckle H, Brock JH, et al. Maternal serum alpha-fetoprotein measurement in antenatal screening for anencephaly and spina bifida in early pregnancy. Report of the UK collaborative study on alpha-fetoprotein in relation to neural-tube defects. *Lancet*. 1977;1(80261):1323–1332.
10. Amniotic-fluid alpha-fetoprotein measurement in antenatal diagnosis of anencephaly and open spina bifida in early pregnancy. Second report of the U.K. Collaborative Study on Alpha-fetoprotein in Relation to Neural-tube Defects. *Lancet*. 1979;2(8144):652–662.
11. Estimating an individual's risk of having a fetus with open spina bifida and the value of repeat alpha-fetoprotein testing. Fourth report of the UK collaborative study on alpha-fetoprotein in relation to neural tube defects. *J Epidemiol Community Health*. 1982;36(2):87–95.

12. Johnson AM, Palomaki GE, Haddow JE. Maternal serum alpha-fetoprotein levels in pregnancies among black and white women with fetal open spina bifida; a United States collaborative study. *Am J Obstet Gynecol*. 1990;162(2):328–331.
13. Brock DJ. The prenatal diagnosis of neural tube defects. *Obstet Gynecol Surv*. 1976;31(1):32–40.
14. Knight GJ. Maternal serum α -fetoprotein screening. In: Hommes FA, ed. *Techniques in Diagnostic Human Biochemical Genetics: A Laboratory Manual*. New York, NY: Wiley-Liss Inc; 1991: 491–518.
15. Burton BK. Elevated maternal serum alpha-fetoprotein (MSAFP): interpretation and follow-up. *Clin Obstet Gynecol*. 1988;31(2):293–305.
16. Palomaki GE, Hill LE, Knight GJ, et al. Second-trimester maternal serum alpha-fetoprotein levels in pregnancies associated with gastroschisis and omphalocele. *Obstet Gynecol*. 1988;71(6, pt 1):906–909.
17. Seppala M, Rapola J, Huttunen NP, et al. Congenital nephrotic syndrome: prenatal diagnosis and genetic counselling by estimation of amniotic-fluid and maternal serum alpha-fetoprotein. *Lancet*. 1976;2(7977):123–125.
18. Abelev GI. Study of the antigenic structure of tumors. *Acta Unio Int Contra Cancrum*. 1963;19:80–82.
19. Tatarinov YS. Detection of embryo-specific alpha-globulin in the blood sera of patients with primary liver cancer. *Vopr Med Khim*. 1964;10:90–91.
20. Abelev GI. Alpha-fetoprotein in ontogenesis and its association with malignant tumors. *Adv Cancer Res*. 1971;14:295–358.
21. Chen DS, Sung JL. Serum alphafetoprotein in hepatocellular carcinoma. *Cancer*. 1977;40(2):779–783.
22. McIntire KR, Vogel CL, Princler GL, Patel IR. Serum alpha-fetoprotein as a biochemical marker for hepatocellular carcinoma. *Cancer Res*. 1972;32(9):1941–1946.
23. Kawai M, Furuhashi Y, Kano T, et al. Alpha-fetoprotein in malignant germ cell tumors of the ovary. *Gynecol Oncol*. 1990;39(2):160–166.
24. Javadpour N. Serum and cellular biologic tumor markers in patients with urologic cancer. *Hum Pathol*. 1979;10(5):557–568.
25. Masopust J, Kithier K, Rádl J, et al. Occurrence of fetoprotein in patients with neoplasms and non-neoplastic diseases. *Int J Cancer*. 1968;3(3):364–373.
26. Waldmann TA, McIntire KR. The use of a radio immunoassay for alpha-fetoprotein in the diagnosis of malignancy. *Cancer*. 1974;34(4 suppl):1510–1515.
27. Javadpour N, McIntire KR, Waldmann TA. Human chorionic gonadotropin (HCG) and alpha-fetoprotein (AFP) in sera and tumor cells of patients with testicular seminoma: a prospective study. *Cancer*. 1978;42(6):2768–2772.
28. Lange PH, McIntire KR, Waldmann TA, et al. Serum alpha fetoprotein and human chorionic gonadotropin in the diagnosis and management of nonseminomatous germ-cell testicular cancer. *N Engl J Med*. 1976;295(22):1237–1240.
29. Javadpour N, McIntire KR, Waldmann TA, et al. The role of the radioimmunoassay of serum alpha-fetoprotein and human chorionic gonadotropin in the intensive chemotherapy and surgery of metastatic testicular tumors. *J Urol*. 1978;119(6):759–762.
30. Kohn J, Orr AH, McElwain TJ, et al. Serum-alpha1-fetoprotein in patients with testicular tumours. *Lancet*. 1976;2(7983):433–436.
31. Perlin E, Engeler JE Jr, Edson M, et al. The value of serial measurement of both human chorionic gonadotropin and alpha-fetoprotein for monitoring germinal cell tumors. *Cancer*. 1976;37(1):215–219.

32. Scardino PT, Cox HD, Waldmann TA, et al. The value of serum tumor markers in the staging and prognosis of germ cell tumors of the testis. *J Urol.* 1977;118(6):994–999.
33. Javadpour N. The role of biologic tumor markers in testicular cancer. *Cancer.* 1980;45(7 suppl):1755–1761.
34. Mason MD. Tumour markers. In: Horwich A, ed. *Testicular Cancer: Investigation and Management.* Baltimore, MD: Williams and Wilkins; 1991:35–50.
35. Toner GC, Geller NL, Tan C, et al. Serum tumor marker half-life during chemotherapy allows early prediction of complete response and survival in nonseminomatous germ cell tumors. *Cancer Res.* 1990;50(18):5904–5910.
36. Kirkpatrick AM, Kirkpatrick KA. Clearance-corrected differencing and other analytic techniques useful in the interpretation of serum AFP values. In: Kirkpatrick AM, Nakamura RM, eds. *Alpha-fetoprotein: Laboratory Procedures and Clinical Applications.* New York, NY: Masson; 1981:135–148.
37. Price P, Hogan SJ, Horwich A. The growth rate of metastatic nonseminomatous germ cell testicular tumours measured by marker production doubling time—I. Theoretical basis and practical application. *Eur J Cancer.* 1990;26(4):450–453.
38. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
39. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
40. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
41. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
42. Seppala M, Ruoslahti E. Radioimmunoassay of maternal serum alpha fetoprotein during pregnancy and delivery. *Am J Obstet Gynecol.* 1972;112(2):208–212.
43. Mukojima T, Hattori N, Nakayama N, et al. Elimination rate of AFP after surgical operation and prognosis of the patients with hepatoblastoma and hepatoma. *Tumor Res.* 1973;8:194–197.
44. Horacek I, Peperell RJ, Hay DL, et al. Detection of fetomaternal haemorrhage by measurement of maternal serum-alpha-fetoprotein. *Lancet.* 1976;2(7978):200.
45. Sizaret P. Equivalence between international units and mass units of alpha-fetoprotein. Report of a collaborative study. *Clin Chimica Acta.* 1979;96(1-2):59–65.
46. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem.* 1999;45(7):942–956.
47. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem.* 1992;38(9):1737–1742.
48. Clinical and Laboratory Standards Institute. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2000. CLSI Document C28-A2.
49. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

50. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
51. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
52. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
53. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP07-ed3.
54. Sizaret P, Anderson SG. The International Reference Preparation for alpha-fetoprotein. *J Biol Stand.* 1976;4(2):149.
55. Sizaret R, Breslow N, Anderson SG. Collaborative study of a preparation of human cord serum for its use as a reference in the assay of alphafoetoprotein. *J Biol Stand.* 1975;3(2):201–223.

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© 2016–2023 Siemens Healthineers.保留所有权利。
Siemens Healthcare Diagnostics Inc.
Siemens Healthcare GmbH