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[siemens.com/healthineers](https://www.siemens.com/healthineers)

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1. Tietz NW, Burlina A, Gerhardt W, et al. Multicenter evaluation of a specific pancreatic isoamylase assay based on a double monoclonal-antibody technique. Clin Chem. 1988;34(10):2096–2102.
2. Cummings ST, Fraser CG. Total amylase and pancreatic isoamylase in serum and urine: Considerations from data on biological variation. Ann Clin Biochem. 1989;26(4):335–340.
3. Gerber M, Naujoks K, Lenz H, Wulff K. A monoclonal antibody that specifically inhibits human salivary alpha-amylase. Clin Chem. 1987;33(7):1158–1162.
4. Svens E, Käpyaho K, Tanner P, Weber TH. Immunocatalytic assay of pancreatic alpha- amylase in serum and urine with a specific monoclonal antibody. Clin Chem. 1989;35(4):662–664.
5. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
6. Clinical and Laboratory Standards Institute. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
7. Clinical and Laboratory Standards Institute. Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
8. Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
9. Siemens Healthcare Diagnostics 存档的数据。
10. Wilding P, Zilva JF, Wilde CE. Transport of specimens for clinical chemistry analysis. Ann Clin Biochem. 1977;14:301-306.
11. Clinical and Laboratory Standards Institute. Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
12. Junge W, Wortmann W, Wilke B, et al. Development and evaluation of assays for the determination of total and pancreatic amylase at 37 ° C according to the principle recommended by the IFCC. Clin Biochem. 2001;34(8):607-615.
13. Junge W, Wortmann W, Wilke B, et al. Erratum: Development and evaluation of assays for the determination of total and pancreatic amylase at 37 ° C according to the principle recommended by the IFCC. Clin Biochem. 2003;36(3):161.
14. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

15. Clinical and Laboratory Standards Institute. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
16. Clinical and Laboratory Standards Institute. Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
17. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.

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Siemens Healthcare Diagnostics Inc.

住所：511 Benedict Avenue, Tarrytown, NY, 10591, United States

生产地址：55 Diamond Road, Crumlin, County Antrim, BT29 4QY,

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© 2019-2020 Siemens Healthcare Diagnostics。保留所有权利。