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当前版本和日期 ^a	Rev. 05, 2023-03	
产品名称	Atellica IM Luteinizing Hormone (LH)	REF 10995635 (110 次测试) REF 10995634 (550 次测试)
产品名称缩写	Atellica IM LH	
测试名称/ID	LH	
系统	Atellica IM Analyzer	
自备材料	Atellica IM CAL B	REF 10995503 (2-包) REF 10995504 (6-包)
可选材料	Atellica IM Multi-Diluent 1	REF 10995637 (2-包) REF 10995638 (6-包) REF 10995639 (小瓶)
	Atellica IM LH MCM	REF 10995636
样本类型	血清、EDTA 血浆、肝素锂血浆	
样本量	25 µL	
测量区间	0.07–200.00 mIU/mL (IU/L)	

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siemens-healthineers.com

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1. Kicklighter EJ, Norman RJ. The gonads. In: Kaplan LA, Pesce AJ, eds. *Clinical Chemistry: Theory, Analysis, Correlation*. 2nd ed. St. Louis, MO: CV Mosby; 1989:650–663.
2. Burtis CA, Ashwood ER. Gonadotropins (follicle-stimulating hormone, luteinizing hormone). In: Tietz NW, ed. *Textbook of Clinical Chemistry*. 2nd ed. Philadelphia: WB Saunders; 1994:1679.
3. Liwnicz BH, Liwnicz RG. The hypothalamopituitary system. In: Kaplan LA, Pesce AJ, eds. *Clinical Chemistry: Theory, Analysis, Correlation*. 2nd ed. St. Louis, MO: CV Mosby; 1989:613–619.
4. Nickel KL. The gonads. In: Kaplan LA, Pesce AJ, eds. *Clinical Chemistry: Theory, Analysis, Correlation*. 3rd ed. St. Louis: CV Mosby, MO; 1996:894.
5. Scott MG, Ladenson JH, Green ED, Gast MJ. Hormonal evaluation of female infertility and reproductive disorders. *Clin Chem*. 1989;35(4):620–629.
6. Butt WR, Blunt SM. The role of the laboratory in the investigation of infertility. *Ann Clin Biochem*. 1988;25(pt 6):601–609.
7. Reissmann T, Felberbaum R, Diedrich K, Engel J, Comaru-Schally AM, Schally AV. Development and applications of luteinizing hormone-releasing hormone antagonists in the treatment of infertility: an overview. *Hum Reprod*. 1995;10(8):1974–1981.
8. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
9. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
10. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
11. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
12. Nerenz RD, Jungheim E, Gronowski AM. Reproductive Endocrinology and Related Disorders. In: Rifai N, Horvath AR, Wittwer C, eds. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 6th ed. Elsevier; 2018:1640-43.
13. Wu A, ed. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. Elsevier; 2006:695.
14. Gooren L. The Neuroendocrine Response of Luteinizing Hormone to Estrogen Administration in Heterosexual, Homosexual, and Transsexual Subjects. *J Clin Endocrinol Metab*. 1986;63(3):583-8.
15. Matsumoto AM. Effects of Chronic Testosterone Administration in Normal Men: Safety and Efficacy of High Dosage Testosterone and Parallel Dose-Dependent Suppression of Luteinizing Hormone, Follicle-Stimulating Hormone, and Sperm Production. *J Clin Endocrinol Metab*. 1990;70(1):282-7.
16. Messinis IE, Templeton A, Baird DT. Endogenous Luteinizing Hormone Surge during Superovulation Induction with Sequential Use of Clomiphene Citrate and Pulsatile Human Menopausal Gonadotropin. *J Clin Endocrinol Metab*. 1985;61(6):1076-80.
17. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45(7):942–956.

18. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem.* 1992;38(9):1737–1742.
19. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
20. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.
21. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
22. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
23. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

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© 2017–2022 Siemens Healthineers。保留所有权利。
Siemens Healthcare Diagnostics Inc.
Siemens Healthcare GmbH