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- Burtis C, Ashwood E. Tietz Fundamentals of Clinical Chemistry, 5th ed. Philadelphia: WB Saunders, 2001:797-799.
- Anderegg G, Flaschka H, Sallmann R, Schwartzenbach G. Metallindikatoren VII ein auf Erdalkalitionen ansprechendes Phtalein und seine analytische Verwendung HELV, Chim Acta 1954; 37:113-120.
- Stern J, Lewis WHP. The colorimetric estimation of calcium in serum with O-cresolphthalein complexone, Clin Chim Acta 1957; 2:576-580.
- Connerty HV, Briggs AR. Determination of serum calcium by means of O-cresolphthalein complexone, Am J Clin Path 1966; 45:290-296.
- Sarker B, Ray C, Chauhan UPS. A new method for determining micro-quantities of calcium in biological materials, Anal Biochem 1967; 20:155-166.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Henry RJ. Clinical Chemistry, Principles and Technics, Harper and Row, New York, NY, 1974, p 656.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1999, p 51 and 1405.
- Henry JB, and LaFever D. An Evaluation of a Calcium Method for the DuPont aca®, DuPont Company, Wilmington, DE (1974).
- Norman PT, Froysa A, Svaland M. Interference of gadodiamide injection (Omniscan®) on the colorimetric determination of serum

calcium. Scand J. Clin Lab Invest 1995; 55; 421-426.

13. Lin J, Idee JM, Port M, Dai A, Berthommier C, Robert M, Raynal I, Devoldere L, Corot C. Interference of magnetic resonance imaging contrast agents with the calcium measurement technique using colorimetric reagents. Journal of Pharmaceutical and Biomedical analysis 21 (1999) 931-943.

14. Barger JD, Dicht JJ. Reference Intervals for Urine on AMY, CA, MG, Na/K, and UP. DuPont Company, Wilmington, DE, October 1983.

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Siemens Healthcare Diagnostics Inc.