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当前版本和日期 ^a	Rev. 02, 2023-05	
产品名称	Atellica CH Apolipoprotein B (APO B)	REF 11097626 (300 次测试)
产品名称缩写	Atellica CH APO B	
测试名称/ID	APO B	
系统	Atellica CI Analyzer	
自备材料	Atellica CH APO A1 & B CAL	REF 11099329
样本类型	血清和血浆 (肝素锂)	
样本量	5 µL	
测量间隔	15–200 mg/dL (0.15–2.00 g/L)	

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siemens-healthineers.com

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1. Alaupovic P. Structure and function of plasma lipoproteins with particular regard to hyperlipoproteinemias and atherosclerosis. *Ann Biol Clin.* 1980;38(2):83–93.
2. Naito HK. The clinical significance of apolipoprotein measurements. *J Clin Immunoassay.* 1986;9:11-20.
3. Burtis CA, Ashwood ER, Bruns DE, eds. *Tietz Fundamentals of Clinical Chemistry.* 6th ed. St. Louis, MO: Saunders; 2008:413–427.
4. Riesen WF, Mordasini R, Salzmann C, Theler A, Gurtner HP. Apoproteins and lipids as discriminators of severity of coronary heart disease. *Atherosclerosis.* 37(1):157–162.
5. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
6. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
7. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
8. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
9. Wu AHB. *Tietz Clinical Guide to Laboratory Tests.* 4th ed. St. Louis, MO: Saunders; 2006:146.
10. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
11. Tietz NW. *Clinical Guide to Laboratory Tests.* 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:68-71.
12. Brown SA, Hutchinson R, Morrisett J, et al. Plasma lipid, lipoprotein cholesterol, and apoprotein distributions in selected US communities. *Arterioscler Thromb.* 1993;13(8):1139-1158.
13. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
14. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
15. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.
16. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
17. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.
18. Siemens Healthcare Diagnostics 存档的数据。

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© 2023 Siemens Healthineers。保留所有权利。
Siemens Healthcare Diagnostics Inc.
Siemens Healthcare GmbH