IMMULITE® 2000

Systems

IGF-I

IMMULITE® 2000 IGF-I

English

Intended Use: For *in vitro* diagnostic use with the IMMULITE® 2000 Systems Analyzers — for the quantitative measurement of insulin-like growth factor I (IGF-I) in serum or heparinized plasma, as an aid in the evaluation of growth disorders.

Catalog Number: **L2KIGF2** (200 tests)

Test Code: **IGF** Color: **Dark Blue**

Summary and Explanation

Insulin-like growth factor I is a single polypeptide chain with three intra-molecule disulfide bonds. It consists of 70 amino acid residues with a molecular mass of 7649 daltons.¹ It is structurally homologous to IGF-II and insulin. IGF-I circulates primarily in a high molecular weight tertiary complex with IGF-binding protein-3 (IGFBP-3) and acid-labile subunit.^{2,3} *In vivo*, IGF-I synthesis is stimulated by growth hormone (GH) and nutritional intake.⁴

In humans, plasma IGF-I levels are barely detectable at birth, rise gradually during childhood, peak during mid-puberty until approximately 40 years of age, then decline gradually.⁵ Maternal plasma levels increase during pregnancy. In the diagnosis of growth disorders, measurements of IGF-I are a useful indicator of growth hormone (GH) secretion. A normal plasma or serum IGF-I concentration is strong evidence against GH deficiency. A low IGF-I value implies GH deficiency and requires additional testing to determine whether GH secretion is subnormal. Measurement of IGF-I is also useful in assessing change of nutritional status.

Measurement of IGF-I in serum is complicated by the presence of acid-labile components and binding proteins.¹ Acid treatment is necessary to release IGF-I to ensure accurate quantitation.^{1,6}

Principle of the Procedure

IMMULITE 2000 IGF-I is a solid-phase, enzyme-labeled chemiluminescent immunometric assay. Sample pretreatment is accomplished in an on-board dilution step.

Incubation Cycles: 1 × 60 minutes

Specimen Collection

Because EDTA would affect results, it should not be used as an anticoagulant.

The use of an ultracentrifuge is recommended to clear lipemic samples.

Hemolyzed or grossly contaminated samples may give erroneous results.

Centrifuging serum samples before a complete clot forms may result in the presence of fibrin. To prevent erroneous results due to the presence of fibrin, ensure that complete clot formation has taken place prior to centrifugation of samples. Some samples, particularly those from patients receiving anticoagulant therapy, may require increased clotting time.

Blood collection tubes from different manufactures may yield differing values, depending on tube materials and additives, including gel or physical barriers, clot activators and/or anticoagulants. IMMULITE 2000 IGF-I has not been tested with all possible variations of tube types. Consult the section on Alternate Sample Types for details on tubes that have been tested.

Volume Required: 20 µL serum or heparinized plasma

Storage: Samples stored at 2–8°C are stable for up to 24 hours. Samples stored at –25°C are stable for 12 months.⁸

Automatic Predilution Factor: 10

Warnings and Precautions

For *in vitro* diagnostic use.

For Professional Use.

H317
P280,
P302+P352,
P333+P313,
P362+P364

Warning!

May cause an allergic skin reaction.

Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention. Take off contaminated clothing and wash it before reuse.

Contains:

2-methylisothiazol-3(2H)-one; IGF-I Reagent Wedges A and B, IGF-I Adjustors, IGF-I Sample Diluent

CAUTION: This device contains material of animal origin and should be handled as a potential carrier and transmitter of disease.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare professional.

Reagents: Store at 2–8°C. Dispose of in accordance with applicable laws.

Safety data sheets (MSDS/SDS) available on siemens-healthineers.com.

Follow universal precautions, and handle all components as if capable of transmitting infectious agents. Source materials derived from human blood were tested and found nonreactive for syphilis; for antibodies to HIV 1 and 2; for hepatitis B surface antigen; and for antibodies to hepatitis C.

Chemiluminescent Substrate: Avoid contamination and exposure to direct sunlight. (See insert.)

Water: Use distilled or deionized water.

Dispose of hazardous or biologically contaminated materials according to the practices of your institution. Discard all materials in a safe and acceptable manner and in compliance with prevailing regulatory requirements.

Materials Supplied

Components are a matched set. Labels on the inside box are needed for the assay.

IGF-I Bead Pack (L2IGF12)

With barcode. 200 beads, coated with monoclonal murine anti-IGF-I. Stable at 2–8°C until expiration date.

L2KIGF2: 1 pack

IGF-I Reagent Wedge (L2IGFA2)

With barcode. 6.5 mL alkaline phosphatase (bovine calf intestine) conjugated to polyclonal rabbit anti-IGF-I in buffer and 6.5 mL IGF-II in buffer. Stable at 2–8°C until expiration date.

L2KIGF2: 1 wedge

Before use, tear off the top of the label at the perforations, without damaging the barcode. Remove the foil seal from the top of wedge; snap the sliding cover down into the ramps on the reagent lid.

IGF-I Adjustors (LGFL, LGFH)

Two vials (Low and High), 4 mL each, of IGF-I in a protein/buffer matrix. (No predilution is required.) Stable at 2–8°C for 30 days after opening or for 6 months (aliquotted) at –20°C.

L2KIGF2: 1 set

Before making an adjustment, place the appropriate Aliquot Labels (supplied with the kit) on test tubes so that the barcodes can be read by the on-board reader.

IGF-I Sample Diluent (L2GFZ)

For the pre-treatment of samples, which is accomplished in the dilution well, through the on-board dilution process. Not intended for diluting samples exceeding the calibration range. 20 mL of concentrated (ready-to-use) IGF-I-free buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months at –20°C.

L2KIGF2: 1 vial

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2KIGF2: 3 labels

Material Required but Not Provided

IGF-I Sample Diluent (L2GFZ)

For the pre-treatment of samples, which is accomplished in the dilution well, through the on-board dilution process. Not intended for diluting samples exceeding the calibration range. 20 mL of concentrated (ready-to-use) IGF-I-free buffer matrix. Stable at 2–8°C for 30 days after opening, or for 6 months at –20°C.

L2GFZ: 1 vial

Barcode labels are provided for use with the diluent. Before use, place an appropriate label on a 16 × 100 mm test tube, so that the barcodes can be read by the on-board reader.

L2GFZ: 3 labels

L2SUBM: Chemiluminescent Substrate

L2PWSM: Probe Wash

L2KPM: Probe Cleaning Kit

LRXT: Reaction Tubes (disposable)

L2ZT: 250 Sample Diluent Test Tubes (16 × 100 mm)

L2ZC: 250 Sample Diluent Tube Caps

LGCOCM: Bi-level IGF control module in a protein/buffer matrix

Also Required

Distilled or deionized water; test tubes; controls

Assay Procedure

Note that for optimal performance, it is important to perform all routine maintenance procedures as defined in the IMMULITE 2000 Systems Operator's Manual.

See the IMMULITE 2000 Systems Operator's Manual for: preparation, setup, dilutions, adjustment, assay and quality control procedures.

Recommended Adjustment Interval:
2 weeks

Quality Control Samples: Use controls or sample pools with at least two levels (low and high) of IGF-I.

Additional quality control material can be used at the discretion of the laboratory. Use the quality control material in accordance with the quality control instructions for use.

In addition, perform quality control:

- Following a valid calibration
- With use of a new lot of reagent
- When troubleshooting test results that do not match clinical conditions or symptoms

Follow government regulations or accreditation requirements for quality control frequency. Individual laboratory quality control programs and procedures may require more frequent quality control testing.

Acceptable performance is achieved when the analyte values obtained are within the expected control interval for the system, as indicated by the manufacturer of the control material or within the interval determined by an internal laboratory quality control procedure.

Follow your laboratory's quality control procedures if the results obtained do not fall within the acceptable limits. For information about entering quality control definitions, refer to the system online help.

Expected Values


A reference range study was performed with the IMMULITE 2000 IGF-I kit on 1804 pediatric and adult samples from an apparently healthy population, yielded the following results, presented in the tables and graphs below.

IGF-I Pediatric Reference Ranges

IGF-I, ng/mL			
Age, y	Gender	Median	Central 95% Range
0–3	Male	44	< 15–129
4–6	Male	96	22–208
7–9	Male	132	40–255
10–11	Male	177	69–316
12–13	Male	305	143–506
14–15	Male	322	177–507
16–18	Male	284	173–414
0–3	Female	68	18–172
4–6	Female	105	35–232
7–9	Female	139	57–277
10–11	Female	248	118–448
12–13	Female	323	170–527
14–15	Female	317	191–496
16–18	Female	291	190–429


IGF-I Pediatric Data

IMMULITE 2000 results for
 $n = 701$ pediatric samples

IGF-I Adult Data

IMMULITE 2000 results for
 $n = 619$ adult samples

IGF-I Adult Reference Ranges

IGF-I, ng/mL

Age, y	Median	Central 95% Range
19–21	207	117–323
22–24	175	99–289
25–29	160	84–259
30–34	136	71–234
35–39	126	63–223
40–44	122	58–219
45–49	120	53–215
50–54	108	48–209
55–59	108	45–210
60–64	112	43–220
65–69	110	40–225
70–79	92	35–216
80–90	94	31–208

IGF-I Tanner Stage Reference Ranges

Female, $n = 235$

IGF-I, ng/mL

Tanner Stage	Median	Central 95% Range
1	186	71–394
2	288	122–508
3	329	164–545
4	319	174–480
5	274	169–400

Male, $n = 249$

IGF-I, ng/mL

Tanner Stage	Median	Central 95% Range
1	144	63–271
2	240	114–411
3	298	166–510
4	290	170–456
5	257	161–384

Consider these limits as *guidelines* only.
Each laboratory should establish its own
reference ranges.

Limitations

Heterophilic antibodies in human serum can react with the immunoglobulins included in the assay components causing interference with *in vitro* immunoassays. [See Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem 1988;34: 27–33.] Samples from patients routinely

exposed to animals or animal serum products can demonstrate this type of interference potentially causing an anomalous result. These reagents have been formulated to minimize the risk of interference; however, potential interactions between rare sera and test components can occur. For diagnostic purposes, the results obtained from this assay should always be used in combination with the clinical examination, patient medical history, and other findings.

The drug asfotase alfa (STRENSIQ®), a recombinant form of alkaline phosphatase, is expected to interfere with *in vitro* diagnostic assays utilizing an alkaline phosphatase detection system. Test samples from patients taking asfotase alfa with a non-alkaline phosphatase methodology.¹⁴

Performance Data

See Tables and Graphs for data *representative* of the assay's performance. Results are expressed in ng/mL. Unless otherwise noted, all were generated on serum samples.

Conversion Factor: ng/mL $\times$ 1 $\rightarrow$ μ g/L
ng/mL $\times$ 0.13 $\rightarrow$ nmol/L

Reportable Range: 15–1000 ng/mL
(WHO NIBSC 1st IS 02/254)

Detection Limits:

Limit of Detection (LoD) = 13.3 ng/mL

Limit of Quantitation (LoQ) = 24.9 ng/mL

LoD and LoQ were determined in accordance with Clinical and Laboratory Standards Institute (CLSI) EP17-A2.¹⁰

High-dose Hook Effect: None up to 90,000 ng/mL

Precision: Samples were processed in duplicate over the course of 20 days, two runs per day, for a total of 40 runs and 80 replicates per CLSI guideline EP05-A3.¹¹ (See "Precision" table.)

Based on internal testing, the overall reproducibility is estimated to be $\leq$ 10% CV for samples tested and includes multiple reagent lots, instruments, days, operators, and replicates. Performance of the assay at individual laboratories may vary.

Linearity: Linearity was assessed in accordance with CLSI EP06-A.¹² The

assay was demonstrated to be linear from 15–1000 ng/mL.

Recovery: Spike recovery studies were carried out using WHO 02/254 spiked into native patient samples. (See "Recovery" table for representative data.)

Specificity: The antibody is highly specific for IGF-I, as assessed per CLSI guideline EP07-A2.¹³ (See "Specificity" table.)

Bilirubin: Presence of conjugated and unconjugated bilirubin in concentrations up to 200 mg/L has no effect on results, within the precision of the assay. Assessed per CLSI guideline EP07-A2.¹³

Hemolysis: Presence of hemoglobin in concentrations up to 500 mg/dL has no effect on results, within the precision of the assay. Assessed per CLSI guideline EP07-A2.¹³

Lipemia: Presence of triglycerides in concentrations up to 3000 mg/dL has no effect on results, within the precision of the assay.

Alternate Sample Types: To assess the effect of alternate sample types, blood was collected from 43 volunteers into plastic serum tubes, plastic heparinized tubes, and plastic gel barrier tubes (SST®). All tubes were from Becton Dickinson. The samples were assayed by the IMMULITE 2000 IGF-I procedure. Results were analyzed by Passing and Bablok linear regression. Because EDTA would affect results, it should not be used as an anticoagulant.

(Heparin Plastic) = 1.04 (Serum Plastic) – 5.3 ng/mL
(SST) = 1.02 (Serum Plastic) – 1.6 ng/mL

Method Comparison 1:

The IMMULITE 2000 IGF-I assay was compared to a commercially available assay for IGF-I (IDS iSYS) on 164 patient samples. (Concentration range: approximately 23–900 ng/mL. See graph.) Results were analyzed by Passing and Bablok linear regression.

(IML 2000) = 0.84 (IDS iSYS) + 1.9 ng/mL

Method Comparison 2:

The IMMULITE 2000 IGF-I assay was compared to the re-standardized IMMULITE 1000 IGF-I on 163 patient samples. (Concentration range: approximately 17–732 ng/mL. See graph.)

Results were analyzed by Passing and Bablok linear regression.

(IML 2000) = 1.01 (IML) + 5.2 ng/mL

References

1) Daughaday WH, Rotwein P. Insulin-like growth factors I and II. Peptide, messenger ribonucleic acid and gene structures, serum, and tissue concentrations. *Endocr Rev* 1989;10: 68-91. 2) Baxter RC, Martin JL, Beniac VA. High molecular weight insulin-like growth factor binding protein complex. Purification and properties of the acid-labile subunit from human serum. *J Biol Chem* 1989;264:11843-8. 3) Rechler MM. Insulin-like growth factor binding proteins. *Vitam Horm* 1993;47:1-114. 4) Lee PD, Rosenfeld RG. Clinical utility of insulin-like growth factor assays. *Pediatrician* 1987;14: 154-61. 5) Hall K, Enberg G, Ritzen M, et al. Somatomedin A levels in serum from healthy children and from children with growth hormone deficiency or delayed puberty. *Acta Endocrinol (Copenh)* 1980;94:155-65. 6) Underwood LE, Murphy MG. Radioimmunoassay of the somatomedins. In: Patrono C, Peskar BA, editors. *Radioimmunoassay in basic and clinical pharmacology. Handbook of experimental pharmacology. vol. 82.* New York: Springer-Verlag; 1987: 561-74. 7) National Committee for Clinical Laboratory Standards. Procedures for the collection of diagnostic blood specimens by venipuncture; approved standard. 4th ed. NCCLS Document H3-A4, Wayne, PA: NCCLS, 1998. 8) Elmlinger MW, Zwirner, M, Kühnel, W. Stability of Insulin-Like Growth Factor (IGF)-I and IGF Binding Protein (IGFBP)-3 Measured by the Immulite® Automated Chemi-Luminescence Assay System in Different Blood Specimens. *Clin Lab* 2005; 51: 145-152. 9) Elmlinger MW, Kühnel W, Weber MM, Ranke MB. Reference ranges for two automated chemiluminescent assays for serum insulin-like growth factor I (IGF-I) and IGF-binding protein 3 (IGFBP-3). *Clin Chem Lab Med* 2004;42:654-64. 10) EP17-A2: Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition, CLSI, Wayne, PA, 2012. 11) EP05-A3: Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition, CLSI, Wayne, PA, 2014. 12) EP06-A: Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline, CLSI, Wayne, PA, 2003 13) EP07-A2: Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition, CLSI, Wayne, PA, 2005. 14) STRENSIQ® (asfotase alfa) injection [package insert]. Boston, MA: Alexion Pharmaceuticals, Inc.; 2020. US License Number: 1743, 2018 (700859).

Technical Assistance

According to EU regulation 2017/746, any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

In the United States, contact Siemens Healthineers Technical Services department. Tel: 877.229.3711. Outside the United States, contact your National Distributor.

siemens-healthineers.com

The Quality System of Siemens Healthcare Diagnostics Products Ltd. is certified to ISO 13485.

Shading indicates technical content that differs from the previous version.

Tables and Graphs

Precision (ng/mL)

	Mean ³	Repeatability ¹		Within-Lab ²	
		SD ⁴	CV ⁵	SD	CV
1	56	3.5	6.3%	4.2	7.6%
2	96	3.2	3.4%	4.6	4.8%
3	200	6.1	3.1%	7.8	3.9%
4	443	13.1	3.0%	16.2	3.7%
5	621	15.8	2.5%	21.3	3.4%
6	863	20.5	2.4%	25.6	3.0%

Recovery (ng/mL)

	Solution ¹	Observed ²	Expected ³	% Bias ⁴
1	—	239	—	—
	A	289	287	1%
	B	470	451	4%
	C	660	641	3%
2	—	216	—	—
	A	261	265	-1%
	B	419	419	0%
	C	631	620	2%


	Solution ¹	Observed ²	Expected ³	% Bias ⁴
3	—	308	—	—
	A	359	350	3%
	B	530	513	3%
	C	717	708	1%
4	—	206	—	—
	A	262	256	2%
	B	431	410	5%
	C	642	608	6%
5	—	74	—	—
	A	135	134	1%
	B	304	279	9%
	C	514	484	6%

Specificity


Compound ¹	Added ²	% Cross-reactivity at 0 ng/mL ³	% Cross-reactivity at 60 ng/mL ⁴
Proinsulin (ng/mL)	155,000	ND	0.00%
Insulin (μIU/mL) μIU/mL × 0.042 → ng/mL	18,000	ND	0.07%
LH (mIU/mL)	70,000	ND	0.04%
TSH (mIU/mL)	310	ND	0.00%
IGF-II (ng/mL)	6000	0.26%	0.27%

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IMMULITE is a trademark of Siemens Healthineers.
All other trademarks and brands are the property of their respective owners.
All rights reserved.