

Dimension® RxL Max® Consumables and Parts List



doclib.siemens-healthineers.com

Consumables List

Part Description	REF	SMN
QuikLYTE® Dilution Check (1 x 50 mL)	S640	10444872
QuikLYTE® Sample Diluent (6 x 500 mL)	S635	10444874
IMT Probe Cleaner (1 x 20 mL)	RD704	10445038
QuikLYTE® Flush Solution (3 x 1000 mL)	S630	10445051
Salt Bridge Solution (1 x 150 mL)	D105	10444873
QuikLYTE® Standard A (3 x 1000 mL)	S620	10444878
QuikLYTE® Standard B (1 x 300 mL)	S625	10444881
Reagent Probe Cleaner (1 x 500 mL)	RD702	10445036
Sample Probe Cleaner (1 x 1000 mL)	RD703	10445035
Purified Water Diluent (4 x 3.5 L)	710615901	10444887
Chemistry Wash (1 x 1.7 L)	RD701	10445052
Empty Flex® reagent cartridge (8 cartridges)	DF99	10445039
Cuvette Cartridge (1)	D828	10445042
Printer Paper (4 rolls)	D829	10445050
Dimension® Reaction Vessels (1000 vessels)	RXV1A	10445044
Sample Cup (1.5 mL, 1000 cups)	DSC4	10445041
Small Sample Container (SSC) (1000 containers)	DSC5	10445040

Parts List

Part Description	REF
Assembly IMT Sample Tubing Kit X0 X1 X2	11241148
RMS Turntable Shuttle Cleaning Tool	11241338
Bottle wrench	10873472
Alignment gauge 0.051 mm (0.002 inches)	11348961
Alignment gauge 0.152 mm (0.006 inches)	11348970
Clot Detect Board Assy	11240731

System Material Kit, SMN 11055691, PN 752705.901 Rev. X

Drawer Number	Qty.	Part Description	Reorder Number
NA	NA	RMS Turntable Shuttle Cleaning Tool	11241338
		NOTE: Available for order, not in kit.	
1	1	Thermal Chamber Filter (Package of 12)	10452017
2	1	Tufoil Compu-Lube Lubricant	10456286
	1	Filter Wheel Dual Sensor	10479324
	1	Reagent Probe Tip and Nut Kit	10455793
	1	Sample Probe Tip and Nut Kit	10457734
	2	Piston Head and L/Washer Nut	10457822
3	1	R1 Probe Tubing	10457690
	1	R2 Probe Tubing	10453640
	1	I-IMT Probe Tubing	10457544
4	1	U-Seal Platen	11241146
	1	Reagent Tray Alignment Gauge	10457556
	1	Sample Probe Alignment Gauge	10457202
5	2	Optical Detector	10452654
	1	Optical Detector (Photometer Home Position ONLY)	10800091
	1	Cuvette Window Extractor Tool	10453606
	1	Cuvette Window (Package of 5)	10455770
6	1	RxL IMT Port/Drain Alignment Gauge	10457561
7	5	Wire Wraps	10451979
	3	E-Ring Push-On	10450182

Drawer Number	Qty.	Part Description	Reorder Number
8	2	Diaphragm	10455823
		NOTE: Customer will receive 3 diaphragms when reordering.	
9	1	Inline Coupling	10451905
	1	Reagent Drain Align Gauge	10457803
10	2	Compression Spring	10452026
	1	Rubber Tip	10453633
	2	E-Ring (3/16 inch)	10452137
11	1	Pump Solenoid Valve	10457244
12	1	Krytox Grease	10455825
	3	Water Filter Adapter	10451992
13	3	Sample Drain Tee	10452108
	1	Sample Drain Brush (Package of 6)	10452209
14	1	Monopump Piston Home Sensor	10452535
	3	E-Ring Push-On	10450183
15	1	Fuse 2.0A (Package of 5)	10457551
	1	Fuse 6.3A (Package of 5)	10457736
	1	Fuse 0.5A (Package of 5)	10457741
	1	Fuse 15.0A (Package of 5)	10457552
	1	Fuse 12.5A (Package of 5)	10457292
	1	Fuse 20A (Package of 5)	10457659
16	1	Fuse 1.6A (Package of 5)	10457738
	1	Fuse 3.1A (Package of 5)	10457316
	1	Cover Thumbscrew	10450379
17	1	Flex [®] Presence Sensor	11241147
18	1	Assembly Kit Sensor Bracket	10457821
19	1	Sample and Aliquot Wheel Lid Optical Sensor	10457190
	1	Syringe, Shoulder Screw (Package of 2)	10457231
20	1	Reflective Switch	10457537
21	.1 Foot	Heat Shrinking Tubing	10450041
	4	Monopump and Rotary Valve Seal	10455826

Drawer Number	Qty.	Part Description	Reorder Number
	1	Reagent Probe Alignment Gauge	10457774
	5	Retaining Ring	10450105
22	1	IMT Diagnostic Cartridge	11240639
	1	IMT Pogo Pins (Package of 8)	10457887
23	1	Assy Flex Tray Sensor	11241301
	1	Waste Manifold Barb Fitting	10805530
	1	Heat Torch	11241145
24	1	Cuvette Wheel Sensor	10456280
	1	Pump Wheel Sensor	10453681
	1	I-IMT/ISE Stylet (Package of 10)	10455798
	1	Felt Replacement Vacuum Muffler	10449794
	2	Tube Fill Guide	10457545
25	1	IMT Pump Tubing Kit (Without TCO2)	10453689
	1	Sample Wheel Alignment Gauge	10457241
	1	Ultrasonics Test Cable	10457710
	1	Alathon Tubing 10 foot roll (0.060 inch I.D.)	10456289
26	1	Sample Tubing Kit	10457692
	1	UPS Bypass Cable	10453656
N/A	1	Aliquot Wheel Alignment Gauge (Located within 750760.901 Canister Package)	10457562

Trademark Information

Dimension, Flex, QuikLYTE, and RxL Max are trademarks of Siemens Healthineers.

All other trademarks and brands are the property of their respective owners.